

Product datasheet for **RG201907**

C9orf6 (FAM206A) (NM_017832) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C9orf6
Synonyms:	C9orf6; CG-8; FAM206A; Simiate
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG201907 representing NM_017832
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTACCGAGCCGAAGCCGCGGAGCCGGTGGTGCCTTCGCTCGTGGATCGATACTTCACTCGCTGGT
ACAAACCGGATGTCAAAGGAAAATTTGTGAGGACCACTGTATACTACAGCACTCTAACCGAATATGTGT
CATCACATTGGCAGAACTCATCCAGTTCTTCAAAGTGAAAAACAATTTAAAGCATTTCTATCAGATC
AGTACCAACTGTAGCAGACTTCAGAACAGGTCTCTGGGAAATTTAAGCGGGGGGCACAGTTTCTAACAG
AGCTTGCACCTCTCTGTAAGATTTACTGCTCAGATGGTGAAGAATATACTGTGTCTAGTTGTGTAGAGG
ACGTTTGATGGAAGTGAATGAAACATCCTCCATAAGCCATCTATTCTTCAAGAAAGCCATCTACTGAA
GGCTACATTGCAGTTGTGTACCCAAATTTGAAGAAAGTAAAGCATAACAGAAGGGTTACTGACACAAA
AACAAATGAAGAAGTCATGGTGAACGCATTATGCCACAACAGCTACGTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG201907 representing NM_017832
Red=Cloning site Green=Tags(s)

MATEPEAAEPVPSLVDRYFTRWYKPDVKGKFCEDHCILQHSNRICVITLAESHPVLQSGKTIKISISYQI
STNCSRLQNKVSGKFKRGAQFLTELAPLCKIYCSDGEEYTVSSCVGRRLMEVNIILHKPSILQEKPSTE
GYIAVVLPKFEESKSITEGLLTQKQYEEVMVKRINATTATS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_017832

ORF Size: 543 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_017832.4](#)

RefSeq Size:	2142 bp
RefSeq ORF:	546 bp
Locus ID:	54942
UniProt ID:	<u>Q9NX38</u>
Cytogenetics:	9q31.3
Gene Summary:	Actin-binding protein that regulates actin polymerization, filopodia dynamics and increases the branching of proximal dendrites of developing neurons.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG201907