

Product datasheet for **RG201909**

LIME (LIME1) (NM_017806) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIME (LIME1) (NM_017806) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	LIME
Synonyms:	dj583P15.4; LIME
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG201909 representing NM_017806 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGCTGCCAGTGTCTGGGCCCTCTGCCCTCTGGGTTCTAGGGTGCTGCGCCCTGCTCCTCTCGC
TGTGGGCGCTGTGCACAGCCTGCCGAGGCCGAGGACGCTGTAGCCCCAGGAAGAGGGCGGGAGGCA
GCGGGCGAGGCTGCAGGGCAGTGCAGCGCGCGGAAGCGTCCCTACTGAGGCGGACCCACCTCTGCTCC
CTCAGCAAGTCGGACACCAGACTGCACGAGCTGCACGGGGCCCGCAGCAGCAGGGCCCTGCGGCCTG
CCAGCATGGATCTCCTGCGCCCACTGGCTGGAGGTGTCCAGGGACATACCGGACCGCAGGCAGCCCC
CTCTGCCTTCCACACCAGGAGCTGCCCGGGCTCTGCCGCGAGCTGCAGCCACCGCAGGGTGCGCTGGC
CTCGAGGCCACCTATTCCAACGTGGGGCTGGCGGCCCTTCCCGGGGTGAGCCTGGCGGCCAGCCCTGTGG
TGGCCGAGTATGCCCAGCTCCAGAAGCGCAAAGGGACCCATCGCAGTCCCCAAGAGCCACAGCAGGGGAA
GACTGAGGTGACCCGCGCGCTCAGGTGGACGTCTGTACTCCAGGGTCTGCAAGCCTAAAGGAGGGAC
CCAGGACCCACACAGACCCGCTGGACCCCAAGGGCCAGGAGCGATTCTGGCCCTGGCGGGTGACCTGG
CCTACCAGACCTCCCGCTCAGGGCCCTGGATGTGGACAGCGGCCCTGGAAAACGTGTATGAGAGCAT
CCGGGAGCTGGGGACCTGCTGGCAGGAGCAGACGTGCGGGGCTGGGACGCCCCCTGCTTCCAGCTGC
CCCAGCCTAGGGAGGGGCTGGAGACCCCTCCCTGCCTCCCTGCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG201909 representing NM_017806
 Red=Cloning site Green=Tags(s)

MGLPVSWAPPALWVLGCCALLLSLWALCTACRRPEDAVAPRKRRRQRARLQGSATAAEASLLRRTHLCS
 LSKSDTRLHELHRGPRSSRALRPASMDLLRPHWLEVS RDITGPQAAPSAFPHQELPRALPAAAATAGCAG
 LEATYSNVGLAALPGVSLAASPVAEYARVQKRKGTHRSPQEPQQGKTEVTPAAQVDVLYSRVCKPKRRD
 PGPTTDPKQGGAILALAGDLAYQTLPLRALDVDSGPLENVYESIRELGDPAGRSSTCGAGTPPASSC
 PSLGRGWRPLPASLP

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_017806

ORF Size: 885 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017806.4</u>
RefSeq Size:	1196 bp
RefSeq ORF:	888 bp
Locus ID:	54923
UniProt ID:	<u>Q9H400</u>
Cytogenetics:	20q13.33
Protein Families:	Transmembrane
Gene Summary:	This gene encodes a transmembrane adaptor protein that links the T and B-cell receptor stimulation to downstream signaling pathways via its association with the Src family kinases Lck and Lyn, respectively. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Mar 2015]

Product images:



Circular map for RG201909