

## Product datasheet for **RG201946**

### ZNRD2 (NM\_006396) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNRD2 (NM\_006396) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ZNRD2  
**Synonyms:** p27; SSSCA1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG201946 representing NM\_006396  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCCTGAACGGAGCTGAAGTCGACGACTTCTCCTGGGAGCCCCGACTGAGGCGGAGACGAAGGTGC  
TGCAGGCGCGACGGGAGCGGCAAGATCGCATCTCCCGGCTCATGGGCGACTATCTGCTGCGCGTTACCG  
CATGCTGGGCGAGACGTGTGCGGACTGCGGGACGATCCTCCTCCAAGACAAACAGCGGAAAATCTACTGC  
GTGGCTTGT CAGGA ACTCGACTCAGACGTGATAAAGATAATCCCCTCTGAATGCCAGGCTGCCCTCT  
CCCAAGCTCGGGAGCACCAGCTGGCCTCAGCCTCAGAGCTCCCCCTGGGCTCTCGACCTGCGCCCCAGCC  
CCCAGTACCTCGTCCGGAGCACTGTGAGGGAGCTGCAGCAGGACTCAAGGCAGCCAGGGGCCACCTGCT  
CCTGCTGTGCCTCCAAATACAGATGTCATGGCCTGCACACAGACAGCCCTCTTGCAAGCTGACCTGGG  
CCTCTGCTGAACTGGGCTCTAGCACCTCCCTGGAGACTAGCATCCAGCTGTGTGGCCTTATCCGCGCATG  
TGCGGAGGCCCTGCGCAGCCTGCAGCAGCTACAGCAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG201946 representing NM\_006396  
 Red=Cloning site Green=Tags(s)

MALNGAEVDDFSWEPPEAETKVLQARRERQDRISRLMGDYLLRGYRMLGETCADCGTILLQDKQRKIYC  
VACQELSDVDKDNPALNAQAALSQAREHQLASASELPLGSRPAPQPPVPRPEHCEGAAAGLKAAQGP  
PAVPPNTDVMACTQTALLQKLTWASAELGSSTSLETSIQLCGLIRACAEALRSLQLQH

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_006396

**ORF Size:** 597 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_006396.3](#)

**RefSeq Size:** 661 bp

**RefSeq ORF:** 600 bp

**Locus ID:** 10534


**UniProt ID:** [O60232](#)

**Cytogenetics:** 11q13.1

**Gene Summary:**

This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren's syndrome. Subcellular localization has not yet been established. [provided by RefSeq, Jul 2008]

**Product images:**



Circular map for RG201946