

Product datasheet for **RG202040**

TXNRD1 (NM_003330) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TXNRD1 (NM_003330) Human Tagged ORF Clone
Symbol:	TXNRD1
Synonyms:	GRIM-12; TR; TR1; TRXR1; TXNR
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG202040 representing NM_003330
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCATGTGAGGACGGTCGGGCCCTGGAAGAACGCTCTCGGAATTGGCCGCGGAAACCGATCTGCCCC
 TTGTGTTTGTGAAACAGAGAAAGATAGGGGCCATGGTCCAACCTTGAAGGCTTATCAGGAGGGCAGACT
 TCAAAAGCTACTAAAAATGAACGGCCCTGAAGATCTTCCCAAGTCTATGACTATGACCTTATCATCATT
 GGAGGTGGCTCAGGAGGTCTGGCAGCTGCTAAGGAGGCAGCCAAATATGGCAAGAAGGTGATGGTCTGG
 ACTTTGCTACTCCACCCTCTTGAAGTAGATGGGGTCTCGGAGGAACATGTGTGAATGTGGGTTGCAT
 ACCTAAAAAAGTATGATCAAGCAGCTTTGTTAGGACAAGCCCTGCAAGACTCTCGAAATTATGGATGG
 AAAGTCGAGGAGACAGTTAAGCATGATTGGGACAGAATGATAGAAGCTGTACAGAATCACATTGGCTCTT
 TGAATTGGGGCTACCGAGTAGCTCTGCGGGAGAAAAAAGTCGTCTATGAGAATGCTTATGGGCAATTTAT
 TGGTCCTCACAGGATTAAGGCAACAAATAATAAGGCAAAAGAAAAATTTATTCAGCAGAGAGATTTCTC
 ATTGCCACTGGTAAAGACCACGTTACTTGGGCATCCCTGGTGACAAAGAATACTGCATCAGCAGTGATG
 ATCTTTTCTCCTTGCCTTACTGCCCGGGTAAGACCCTGGTTGTTGGAGCATCCTATGTCGCTTTGGAGTG
 CGCTGGATTTCTTGTGTTATTGGTTTAGACGCTACTGTTATGGTTAGGTCCATTCTTCTTAGAGGATTT
 GACCAGGACATGGCCAACAAAATGGTGAACACATGGAAGAACATGGCATCAAGTTTATAAGACAGTTTCG
 TACCAATTAAGTTGAACAAATGAAGCAGGGACACCAGGCCGACTCAGAGTAGTAGCTCAGTCCACCAA
 TAGTGAGGAAATCATTGAAGGAGAATAATAACGGTGTGCTGGCAATAGGAAGAGATGCTTGCACAAGA
 AAAATTGGCTTAGAAACCGTAGGGTGAAGATAAATGAAAAGACTGAAAAAATACCTGTCACAGATGAAG
 AACAGCAATGTGCCTTACATCTATGCCATTGGCGATATATTGGAGGATAAGGTGGAGCTCACCCAGT
 TGCAATCCAGGCAGGAAGATTGCTGGCTCAGAGGCTCTATGCAGGTTCCACTGTCAAGTGTGACTATGAA
 AATGTTCCAACCACTGTATTTACTCCTTTGGAATATGGTCTTGTGGCCTTTCTGAGGAGAAAGCTGTGG
 AGAAGTTTGGGAAGAAAATATTGAGGTTTACCATAGTTACTTTTGGCCATTGGAATGGACGATTCGCTC
 AAGAGATAACAACAAATGTTATGCAAAAATAATCTGTAATACTAAAGACAATGAACGTGTTGTGGCTTT
 CACGTAAGGTTCCAAATGCTGGAGAAGTTACACAAGGCTTTGCAGCTGCGCTCAAATGTGGACTGACCA
 AAAAGCAGCTGGACAGCACAATTGGAATCCACCCTGTCTGTGCAGAGGATTACACAACATTGTCTGTGAC
 CAAGCGCTCTGGGCAAGCATCCTCCAGGCTGGCTGCTGAGGT

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>RG202040 representing NM_003330
 Red=Cloning site Green=Tags(s)

MSCEDGRALEGLSELAETDLPVVFVKQRKIGGHGPTLKAYQEGRLQKLLKMNGPEDLPKSYDYDLII
 GGGSGGLAAAKEAAQYGKKVMVLDVFTPTPLGTRWGLGGTCVNVGCIIPKLMHQALLGQALQDSRNYGW
 KVEETVKHDWDRMIEAVQNHIGSLNWGYRVALREKKVVYENAYGQFIGPHRIKATNNGKKEKIYSAERFL
 IATGERPRYLGIIPGDKEYCISDDLFSLPYCPGKTLVVGASYVALECAFLAGIIGLDVTVMVRISILLRFG
 DQDMANKIGEHEEHGKIFIRQFVPIKVEQIEAGTPGRLRVVAQSTNSEEIIEGEYNTVMLAIGRDACTR
 KIGLETVGVKINEKTGKIPVTDEEQTNVPYIYIAGDILEDKVELTPVAIQAGRLLAQRLYAGSTVKCDYE
 NVPTTVFTPLEYGACGLSEEKAVEKFGREENIEVYHSYFWPLEWTIPSRDNNKCYAKIICNTKDNERVVGF
 HVLGPNAGEVTQGFAAALKCGLTKKQLDSTIGIHPVCAEVFTTLSVTKRSGASILQAGC*G

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_003330

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#) The expression of this clone is not guaranteed due to the nature of selenoproteins.

OTI Annotation: This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003330.4](#)

RefSeq Size: 4206 bp


RefSeq ORF: 1656 bp

Locus ID: 7296

UniProt ID: [Q16881](#)

Cytogenetics:	12q23.3
Domains:	pyr_redox, pyr_redox_dim
Protein Families:	Druggable Genome
Protein Pathways:	Pyrimidine metabolism
Gene Summary:	<p>The protein encoded by this gene belongs to the pyridine nucleotide-disulfide oxidoreductase family, and is a member of the thioredoxin (Trx) system. Three thioredoxin reductase (TrxR) isozymes are found in mammals. TrxRs are selenocysteine-containing flavoenzymes, which reduce thioredoxins, as well as other substrates, and play a key role in redox homeostasis. This gene encodes an ubiquitously expressed, cytosolic form of TrxR, which functions as a homodimer containing FAD, and selenocysteine (Sec) at the active site. Sec is encoded by UGA codon that normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, the Sec insertion sequence (SECIS) element, which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. Alternative splicing, primarily at the 5' end, results in transcript variants encoding same or different isoforms, including a glutaredoxin-containing isoform that is predominantly expressed in testis. [provided by RefSeq, May 2017]</p>

Product images:



Circular map for RG202040