

## Product datasheet for **RG202075**

### **IL8 (CXCL8) (NM\_000584) Human Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	IL8 (CXCL8) (NM_000584) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	IL8
Synonyms:	GCP-1; GCP1; IL8; LECT; LUCT; LYNAP; MDNCF; MONAP; NAF; NAP-1; NAP1; SCYB8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202075 representing NM_000584 Red=Cloning site Blue=ORF Green=Tags(s)

```
TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCCGATCGCC
```

```
ATGACTTCCAAGCTGGCCGTGGCTCTTGGCAGCCTTCCTGATTTCTGCAGCTCTGTGTGAAGGTGCAG  
TTTTGCCAAGGAGTGCTAAAGAACTTAGATGTCAGTGCATAAAGACATACTCCAAACCTTTCCACCCCAA  
ATTTATCAAAGAAGTGAAGTGATTGAGAGTGGACCACACTGCGCCAACACAGAAATTATTGTAAGCTT  
TCTGATGGAAGAGAGCTCTGTCTGGACCCCAAGAAAAGTGGTGCAGAGGGTTGTGGAGAAGTTTTTGA  
AGAGGGCTGAGAATTCA
```

```
ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA
```

Protein Sequence:	>RG202075 representing NM_000584 Red=Cloning site Green=Tags(s)
-------------------	--

```
MTSKLAVALLAAFLISAALCEGAVLPRSAKELRCQCIKTYSKPFHPKFIKELRVIESGPHCANTEIIVKL  
SDGRELCLDPKENWVQRVVEKFLKRAENS
```

```
TRTRPLE - GFP Tag - V
```

Restriction Sites:	Sgfl-MluI
--------------------	-----------



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_000584

**ORF Size:** 297 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**RefSeq:** [NM\\_000584.4](#)

**RefSeq Size:** 1666 bp

**RefSeq ORF:** 300 bp

<b>Locus ID:</b>	3576
<b>UniProt ID:</b>	<a href="#">P10145</a>
<b>Cytogenetics:</b>	4q13.3
<b>Domains:</b>	IL8
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Bladder cancer, Chemokine signaling pathway, Cytokine-cytokine receptor interaction, Epithelial cell signaling in Helicobacter pylori infection, NOD-like receptor signaling pathway, Pathways in cancer, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of the CXC chemokine family and is a major mediator of the inflammatory response. The encoded protein is commonly referred to as interleukin-8 (IL-8). IL-8 is secreted by mononuclear macrophages, neutrophils, eosinophils, T lymphocytes, epithelial cells, and fibroblasts. It functions as a chemotactic factor by guiding the neutrophils to the site of infection. Bacterial and viral products rapidly induce IL-8 expression. IL-8 also participates with other cytokines in the proinflammatory signaling cascade and plays a role in systemic inflammatory response syndrome (SIRS). This gene is believed to play a role in the pathogenesis of the lower respiratory tract infection bronchiolitis, a common respiratory tract disease caused by the respiratory syncytial virus (RSV). The overproduction of this proinflammatory protein is thought to cause the lung inflammation associated with cystic fibrosis. This proinflammatory protein is also suspected of playing a role in coronary artery disease and endothelial dysfunction. This protein is also secreted by tumor cells and promotes tumor migration, invasion, angiogenesis and metastasis. This chemokine is also a potent angiogenic factor. The binding of IL-8 to one of its receptors (IL-8RB/CXCR2) increases the permeability of blood vessels and increasing levels of IL-8 are positively correlated with increased severity of multiple disease outcomes (eg, sepsis). This gene and other members of the CXC chemokine gene family form a gene cluster in a region of chromosome 4q. [provided by RefSeq, May 2020]</p>

Product images:



Circular map for RG202075