

Product datasheet for **RG202141**

ARF1 (NM_001658) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ARF1 (NM_001658) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ARF1
Synonyms: PVNH8
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG202141 representing NM_001658
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGAACATCTTCGCCAACCTCTCAAGGGCCTTTTTGGCAAAAAGAAATGCGCATCCTCATGGTGG
GCCTGGATGCTGCAGGGAAGACCACGATCCTCTACAAGCTTAAGCTGGGTGAGATCGTGACCACCATTC
CACCATAGGCTTCAACGTGGAACCGTGGAGTACAAGAACATCAGCTTCACTGTGTGGGACGTGGGTGGC
CAGGACAAGATCCGGCCCTGTGGCCCACTACTCCAGAACACACAAGGCCTGATCTTCGTGGTGGACA
GCAATGACAGAGAGCGTGTGAACGAGGCCCGTGAAGGAGCTCATGAGGATGCTGGCCGAGGACGAGCTCCG
GGATGCTGTCTCCTGGTTCGCCAACAGCAGGACCTCCCAACGCCATGAATGCGGCCGAGATCACA
GACAAGCTGGGGCTGCACTCACTACGCCACAGGAAGTGGTACATTCAGGCCACCTGCGCCACCAGCGGCG
ACGGGCTCTATGAAGGACTGGACTGGCTGTCCAATCAGCTCCGGAACCAGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG202141 representing NM_001658
Red=Cloning site Green=Tags(s)

MGNIFANLFGKLFGKEMRILMVGLDAAGKTTILYKLLKGEIVTTIPTIGFNVETVEYKNISFTVWDVGG
QDKIRPLWRHYFQNTQGLIFVVDSDNRERVNEAREELMRMLAEDELRDVLLVFANKQDLPNAMAAEIT
DKLGLHSLRHRNWIQATCATSGDGLYEGLDWLSNQLRNQK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:


ACCN: NM_001658

ORF Size: 543 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001658.4](#)

RefSeq Size: 1901 bp

RefSeq ORF: 546 bp


Locus ID: 375

UniProt ID: [P84077](#)

Cytogenetics: 1q42.13
Domains: RAB, SAR, ARF, arf
Protein Pathways: Vibrio cholerae infection

Gene Summary: ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. The gene products, including 6 ARF proteins and 11 ARF-like proteins, constitute a family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2 and ARF3), class II (ARF4 and ARF5) and class III (ARF6), and members of each class share a common gene organization. The ARF1 protein is localized to the Golgi apparatus and has a central role in intra-Golgi transport. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG202141