

Product datasheet for **RG202293**

MFSD3 (NM_138431) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	MFSD3
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG202293 representing NM_138431
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCGCGGGAAGTTGCTGCCGCTGGCCGGCCTCTACCTGGTGCAGGGCCTGCCCTACGGGCTCCAGTCCG
GCCTCTGCCAGTGCTGCTGCGTCCCGCGGCCTCTCGCTGACGCGCTGGGGCTGGCCAAGTTCTGTA
CGCTCCGTGGCTGCTCAAGCTGGCTTGGGCCCCGCTGGTGGACGCGCAGGGCTCGGCCAGGGCCTGGGTG
ACGCGCAGCACGGCGGGCCTGGGCCTGGTGTGTGGGCTGCTTGCCGGGCTGCCCCCTCTGGAGCTGGCC
AGGCCGGGCTGCCCGCCGCTGTGGCGGGTTGCTGCTGTTGTTGAACCTGGGTGCCGCCATGCAGGATGT
GGCCCTGGACGCGCTGGCTGTGCAGCTGCTGGAGCCGGCCGAACCTGGGGCCGGCAATACCGTGCAGGTG
GTCGCGTACAAGCTGGGGGCCGCGCTAGCTGGGGCGCGCTGCTGGCGCTGCTGCCACCTTCTCGTGCG
CGCAACTCTTTCTGCTCCTGGCTGCCACCTACTGGCTGGCCGCGGCCCTGGCCTGGGCTGCACGACCCCT
GCGGCGGCTCCACAGCAGCCCCCTCCGAGCAGCGTCCCCACCCGCGCACCTTCTGCGGGACGTGCTA
GCCGTGCCGGGACCGTGTGGACGGCAGGCTTTGTGCTCACCTACAAGCTGGGTGAGCAGGGTGCCAGCA
GCCTGTTTCTCTTCTCCTGCTGGACCACGCGTTTCTGCTCCCGAGTTGGGACTGTGGAATGGTGTGGG
TGCTGTGGTCTGCTCCATCGCTGGCTCCTCCCTGGGTGGGACCTTGCTGGCCAAGCACTGGAAACTGCTG
CCTCTGTTGAGGTGCGTGTGCGCTTCGCGCTCGGGGGCCTAGCCTGTGAGACTGCCTTGGTCTTCCACC
TGGACACCTTGGGGGCCAGCATGGACGCTGGCACAATCTTGAGAGGGTCAGCCTTGTGAGCCTATGCTC
GCAGCACTTCTTGGGAGGCTGTGTCACACAGTCACCTTCACTGGGATGATGCGCTGCAGCCAGCTGGCC
CCCAGGGCCCTGCAGGCCACACACTACAGCCTTCTGGCCACGCTGGAGCTGCTGGGGAAGCTGCTGCTGG
GCACCTGCGCGGAGGCTGGTGTGTTGGGGCCACATCCCTGCTTCTGCTCCTGCTCATCTCTC
TGCTTCCCGTTCTGTACCTGGACCTAGCACCCAGCACCTTCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG202293 representing NM_138431
Red=Cloning site Green=Tags(s)

MRGKLLPLAGLYLVQGLPYGLQSGLLPVLLRAGGLSLTRVGLAKVLYAPWLLKLAWAPLVDAQGSARAWV
TRSTAGLGLVCGLLAGLPPPGAGQAGLPAAVAGLLLLLNLGAAMQDVALDALAVQLLEPAELGPGNTVQV
VAYKLGAAALAGGALLALLPTFSWPQLFLLLAATYWLAAALAWAAPALRRLPQQPPSEQRPHTAHLRLDVL
AVPGTVWTAGFVLTYYLGEQGASSLFPLLLLDHGVSAPELGLWNGVGAVVCSIAGSSLGGLLAKHWKLL
PLLRSVLRFRLGGLACQTAIVFHLDTLGASMDAGTILRGSALLSLCLQHFLGGLVTTVTFTGMMRCSQLA
PRALQATHYSLLATLELLGKLLGLTAGGLADGLGPHPCFLLLLILSAFPVLYLDLAPSTFL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_138431


ORF Size: 1236 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min.
	2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
	3. Close the tube and incubate for 10 minutes at room temperature.
	4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
	5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_138431.1, NP_612440.1</u>
RefSeq Size:	1468 bp
RefSeq ORF:	1239 bp
Locus ID:	113655
UniProt ID:	<u>Q96ES6</u>
Cytogenetics:	8q24.3
Protein Families:	Transmembrane

Product images:


Circular map for RG202293