

Product datasheet for **RG202534**

TMEM204 (NM_024600) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	TMEM204
Synonyms:	C16orf30; CLP24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG202534 representing NM_024600
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGACCGTGCAGAGACTCGTGGCCGCGGCCGTGCTGGTGGCCCTGGTCTCACTCATCCTCAACAACGTGG
CGGCCTTCACCTCCAACCTGGGTGTGCCAGACGCTGGAGGATGGGCGCAGGCGCAGCGTGGGGCTGTGGAG
GTCTGCTGGCTGGTGGACAGGACCCGGGAGGGCCGAGCCCTGGGGCCAGAGCCGGCCAGGTGGACGCA
CATGACTGTGAGGCGCTGGGCTGGGGCTCCGAGGCAGCCGGCTTCAGGAGTCCCGAGGCACCGTCAAAC
TGCAGTTCGACATGATGCGCGCTGCAACCTGGTGGCCACGGCCGCGCTACCGCAGGCCAGCTCACCTT
CCTCCTGGGGCTGGTGGCCCTGCCCCGTGTCACCCGACGCCCCGTGCTGGGAGGAGGCCATGGCCGCT
GCATTCCAACCTGGCGAGTTTGTCTGGTCATCGGGCTCGTGACTTTCTACAGAATTGGCCCATACACCA
ACCTGTCCTGGTCTGTACCTGAACATTGGCGCTGCCTTCTGGCCACGCTGGCGGCAGCCATGCTCAT
CTGGAACATTCTCCACAAGAGGGAGGACTGCATGGCCCCCGGGTGATTGTCATCAGCCGCTCCCTGACA
GCGCGCTTTCGCCGTGGGCTGGACAATGACTACGTGGAGTCAACATGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG202534 representing NM_024600
Red=Cloning site Green=Tags(s)


MTVQRLVAAAVLVALVSLILNNVAFTSNWVCQTLEDGRRRSVGLWRSCWLVDRTRGPSPGARAGQVDA
HDCEALGWGSEAAGFQESRGTVKLQFDMMRACNLVATAALTAGQLTFLLGLVGLPLLSPDAPCWEEMAA
AFQLASFVLVIGLVTFYRIGPYTNLSWSCYLNIGACLLATLAAAMLIWNILHKREDCMAPRVIVISRSLT
ARFRRLDNDYVESPC

TRTRPLE - GFP Tag - V



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024600

ORF Size: 678 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024600.6](#)

RefSeq Size: 1871 bp

RefSeq ORF: 681 bp

Locus ID: 79652


UniProt ID: [Q9BSN7](#)

Cytogenetics: 16p13.3

Protein Families: Transmembrane

Gene Summary: C16ORF30 plays a role in cell adhesion and cellular permeability at adherens junctions (Kearsey et al., 2004 [PubMed 15206924]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RG202534