

Product datasheet for **RG202554**

C14ORF140 (ZC2HC1C) (NM_001042430) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C14ORF140 (ZC2HC1C) (NM_001042430) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C14ORF140
Synonyms:	C14orf140; FAM164C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202554 representing NM_001042430 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGGTCTCCAGCGGTTGGCGTCACATCTGCCTGTGGCGTTATGCTCCACATAATACAACGGAAG
CTCCAGGGCCCCACTCAGCCAAGCAAGACTCTTACGAACAAGGTGACTCTCCAGCAGTCCTTGAAGGG
GCACCTGAGGAACAATTTCCAGAAGCAGCTTTTGAGCAACAAAGATTGATTCTGGATAAAGTCTATACT
CACCCCAATGGAACACCCAAACAAAAGCCGGAGCTACTCTATCCCCACTGTACTGGAATCAGCCAGC
AAGATCCAGAAAGTGATTCCCAGGGCCAAGGAAATGGTTTGTCTTACTCGTCAGGCCCTCAATCCTGGTA
TCCCAAAGCCAATAACCAGGACTTTATCCCTTTACAAAGAAACGAGTTGGAGTGGACCGGGCGTCCCA
TTGAAACCCATGGTCCACAGGAAGTCGTGACGTACAGGTGAGGCTGGCACTGATGGGACCATAATGTCT
ACCCAAGGCCCTGAGCCGAGAGATTTTCATCTAGGAACTTTGGTGTGAGGAACCAGGGCAACTTTTC
TGTGGTTGGTACTGTTCTTGCTGCCACGCAGGCGGAGAAGGCCGTGGCAAACCTTTGACAGGACGGAGTGG
GTGCAGATCCGAAGACTAGAAGCTGCAGGGGAGAGCTTAGAGGAGGAAATCCGAAGAAAGCAGATTCTCC
TGAGGGGAAAGCTGAAGAAGACAGAGGAGGAACCTCAGAAGGATCCAGACGCAAAAGGAACAGGCCAAGGA
AAATGAAAACGGAGACTACAGAAAATTATACTCCCAGGAGCAGAGTTAAAGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG202554 representing NM_001042430
 Red=Cloning site Green=Tags(s)

MAGLQRLASHLPVGVMLPHNTTEAPGPHSAKQDSYEQGDSSQQSLKGHLRNNFQKQLLSNKEILDKVYT
 HPKWNTQTKARSYSYPHCTGISQQDPESDSQGQGNLFYSSGPQSWYPKANNQDFIPFTKKRVGVDRAFP
 LKPMVHRKSCSTGEAGTDGDHNVYPRPPEPREFSSRNFGVRNQGNSVVGTVLAATQAEKAVANFDRTEW
 VQIRRLAAGESLEEEIRRKQILLRGKLLKTEELRRIQTQKEQAKENENGELQKIILPRSRVKG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001042430

ORF Size: 825 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042430.3](#)

RefSeq Size: 1578 bp

RefSeq ORF: 828 bp

Locus ID: 79696

UniProt ID: [Q53FD0](#)

Cytogenetics: 14q24.3

Product images:



Circular map for RG202554