

Product datasheet for **RG202640**

MAD2L1 binding protein (MAD2L1BP) (NM_014628) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAD2L1 binding protein (MAD2L1BP) (NM_014628) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MAD2L1 binding protein
Synonyms:	CMT2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202640 representing NM_014628 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCGCCGGAGGCGGAGGTTCTGTCTCAGCCGAGTCCCTGATTTGGAGTGGTATGAGAAGTCCG
AAGAACTCACGCCTCCCAGATAGAAGTCTTGTAGACAAGCTCTACGCAGGAACCTCTCAACGCTTCGGA
GGCCTTTGCCAAGAGACTGCATGGTACCAGTGGTGTTCCTGGGCCTGTGAGCCAGGAAGGCTGCTGT
CAGTTTACTTGTGAAGTCTAAAGCATATCATGTATCAACGCCAGCAGCTCCCTCTGCCCTATGAACAGC
TTAAGCACTTTTACCGAAAACCTTCTCCCAGGCAGAGGAGATGCTGAAGAAGAAACCTCGGGCCACCAC
TGAGGTGAGCAGCAGGAAATGCCAACAAGCCCTGGCAGAACTGGAGAGTGTCTCAGCCACCTGGAGGAC
TTCTTTGCACGGACACTAGTACCGCAGTGTCTGATTCTCCTGGGGCAATGCCCTAAGCCCAAGGAGT
TCTATGAAGTTCGACTTGTCTCTGCTGGCCCCCTACAGCGTGGACCAGAGCCTGAGCACAGCAGCTTGTTT
GCGCCGTCTCTCCGAGCCATATTCATGGCTGATGCCTTTAGCGAGCTTCAGGCTCCTCCACTCATGGGC
ACCGTCGTCATGGCACAGGGACCCGAACTGTGGAGAAGATTGGTTTCGACCCAAGCTCAACTATCGAG
TGCCAGCCGGGGCCATAAACTGACTGTGACCCGTGTCATGTGGCAGACCTCCATCCGAACCACGGCTTG
GGAAGACTACATTTGGTTCCAGGCACCAAGTACATTTAAAGGCTTCCGCGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG202640 representing NM_014628
 Red=Cloning site Green=Tags(s)

MAAPEAEVLSSAAVPDLEWYEKSEETHASQIELLETSSSTQEPLNASEAF CPRDCMVPVFPVPGPVSQEGCC
 QFTCELLKHIMYQRQQLPLPYEQLKHFYRKPSQAEEMLKKKPRATTEVSSRKCCQALAELESVLSHLED
 FFARTLVPRVLILLGGNALSPKEFYELDL SLLAPYSVDQSLSTAACLRRLFRAIFMADAFSELQAPPLMG
 TVVMAQGHRCGEDWFRPKLNRYRVP SRGHKLT VTLSCGRPSIRTTAWEDIWFQAPVTFKGFRE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_014628

ORF Size: 822 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014628.3](#)

RefSeq Size: 1283 bp

RefSeq ORF: 825 bp

Locus ID: 9587


UniProt ID: [Q15013](#)

Cytogenetics: 6p21.1

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene was identified as a binding protein of the MAD2 mitotic arrest deficient-like 1 (MAD2/MAD2L1). MAD2 is a key component of the spindle checkpoint that delays the onset of anaphase until all the kinetochores are attached to the spindle. This protein may interact with the spindle checkpoint and coordinate cell cycle events in late mitosis. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG202640