

Product datasheet for **RG202665**

UFM1 (NM_016617) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: UFM1 (NM_016617) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: UFM1
Synonyms: BM-002; C13orf20; HLD14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG202665 representing NM_016617
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTGGAAGTTTCCTTTAAGATCACGCTGACGTCGGACCCACGGCTGCCGTACAAAGTACTCAGTGTTT
 CTGAAAGTACACCTTTCACAGCAGTCTTAAAGTTTGCAGCAGAAGAATTTAAAGTTCTGCTGCAACAAG
 TGCAATTATTACCAATGATGGAATAGGAATAAATCCTGCACAGACTGCTGGAAATGTTTTCTAAACAT
 GGTTCAGAACTGCGGATTATTCCTAGAGATCGTGTGGAAGTTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG202665 representing NM_016617
 Red=Cloning site Green=Tags(s)
 MSKVSFKITLTSDPRLPYKVLSPPESTPFTAVLKFAAEFKVPAATSAIITNDGIGINPAQTAGNVFLKH
 GSELRIIPDRVGSC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_016617

ORF Size: 255 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_016617.4](#)

RefSeq Size: 2529 bp

RefSeq ORF: 258 bp

Locus ID: 51569


UniProt ID: [P61960](#)

Cytogenetics: 13q13.3

Domains: UPF0185

Gene Summary: UFM1 is a ubiquitin-like protein that is conjugated to target proteins by E1-like activating enzyme UBA5 (UBE1DC1; MIM 610552) and E2-like conjugating enzyme UFC1 (MIM 610554) in a manner analogous to ubiquitylation (see UBE2M; MIM 603173) (Komatsu et al., 2004 [PubMed 15071506]).[supplied by OMIM, Dec 2008]

Product images:



Circular map for RG202665