

Product datasheet for **RG202696**

CRSP9 (MED7) (NM_004270) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CRSP9 (MED7) (NM_004270) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CRSP9
Synonyms:	ARC34; CRSP9; CRSP33
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202696 representing NM_004270 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGAACCACAGCAAGTGAGTGCACTTCCACCACCTCCAATGCAATATATCAAGGAATATACGGATG
AAAATATTCAAGAAGGCTTAGCTCCAAGCCTCCCCTCCAATAAAAGACAGTTACATGATGTTGGCAA
TCAGTTCCAATGTGATGATCTTATCATCCGCCCTTGGAAAGTCAGGGCATCGAACGGCTTCATCTATG
CAGTTTGATCACAAGAAAGAAGTGAAGAACTTAATATGTCTATCCTTATTAATTTCTTGGACCTTTT
ATATTTAATAAGGAGCCCTGGGAGTATAAAACGAGAAGAGAACTAGAAGATCTTAAGCTGCTTTTGT
ACACGTGCATCATCTTATAAATGAATACCGACCCACCAAGCAAGAGAGACCTTGAGAGTCATGATGGAG
GTCCAGAAACGTCAACGGCTTGAACAGCTGAGAGATTTCAAAGCACCTGGAACGAGTAATTGAAATGA
TTCAGAATTGCTTGGCTTCTTGCCTGATGATTTGCCTCATTGAGAAGCAGGAATGAGAGTAAAACTGA
ACCAATGGATGCTGATGATAGCAACAATTGACTGGACAGAATGAACATCAAAGAGAAAATTCAGGTCAT
AGGAGAGATCAGATTATAGAGAAAGATGCTGCCTTGTGTGCCTAATTGATGAGATGAATGAAAGACCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG202696 representing NM_004270
 Red=Cloning site Green=Tags(s)

MGEPQQVSALPPPPMQYIKEYTDENIQEGLAPKPPPIKDSYMMFGNQFQCDDLIIIRPLESQGIERLHPM
 QFDHKELRKLNMSILINFLDLLDILIRSPGSIKREKLEDLKLFFVHVHHLINERYPHQARETLRVMMME
 VQKRQRLETAERFQKHLERVIEMIQNCLASLPDDLPHSEAGMRVKTEPMDADDSNNCTGQNEHQRENSGH
 RRDQIIIEKDAALCVLIDEMNERP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_004270

ORF Size: 699 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_004270.5](#)

RefSeq Size: 1066 bp

RefSeq ORF: 702 bp

Locus ID: 9443


UniProt ID: [O43513](#)

Cytogenetics: 5q33.3

Protein Families: Druggable Genome, Transcription Factors

Gene Summary: The activation of gene transcription is a multistep process that is triggered by factors that recognize transcriptional enhancer sites in DNA. These factors work with co-activators to direct transcriptional initiation by the RNA polymerase II apparatus. The protein encoded by this gene is a subunit of the CRSP (cofactor required for SP1 activation) complex, which, along with TFIID, is required for efficient activation by SP1. This protein is also a component of other multisubunit complexes e.g. thyroid hormone receptor-(TR-) associated proteins which interact with TR and facilitate TR function on DNA templates in conjunction with initiation factors and cofactors. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG202696