

## Product datasheet for **RG202724**

### ADTRP (NM\_032744) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ADTRP (NM_032744) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	ADTRP
Synonyms:	AIG1L; C6orf105; dj413H6.1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202724 representing NM_032744 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACGAAGACTTCTACATGCATATACCACTTCCTTGTCTGAGCTGGTATACTTCTCAATTATTACA  
TCTCACAGGAAGGAAAAGACGAGGTGAAACCCAAAATCTTGGCAAATGGTCAAGGTGAAATATATGAC  
GCTGCTTAATCTGCTCTTGCAGACCATTTCTACGGGGTCACCTGCCTGGATGATGTGCTGAAAAGAACC  
AAAGGGGAAAAGACATTAAGTTCCTAACTGCCTCAGAGACCTGCTTTTACCACCTCTGGCTTTTCTCTG  
TATCCACGTTTGTATTTTGGCATTCTGGATCCTTTTCTCTACAATCGAGATCTCATTTACCCCAAGGT  
CCTAGATACTGTCAATCCCGTGTGGCTGAATCATGCAATGCACACTTTCAATTTCCCATCACATTGGCT  
GAAGTCGTCCTCAGGCCTCACTCCTATCCATCAAAGAAGACAGGACTCACCTTGCTGGCTGCTGCCAGCA  
TTGCTTACATCAGCCGCATCCTATGGCTCTACTTTGAGACGGGTACCTGGGTGTATCCTGTGTTTGCCAA  
ACTCAGCCTCTTGGGTCTAGCAGCTTTCTCTCTCAGCTACGTCTTCATCGCCAGCATCTACCTACTT  
GGAGAGAAGCTCAACCACTGAAATGGGGTGACATGAGGCAGCCACGGAAGAAGAGGAAG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG202724 representing NM\_032744  
 Red=Cloning site Green=Tags(s)

MTKTSTCIYHFLVLSWYTFLNYYISQEGKDEVKPKILANGARWKYMTLLNLLQTIFYGVTCLDDVLKRT  
 KGGKDIKFLTAFRDLLFTLAFPVSTFVFLAFWILFLYNRDLIYPKVLDTVIPVWLNHAMHTFIFPITLA  
 EVVLRPHSYPSSKKTGLTLAAASIAIYISRILWL YFETGTWVYPVFAKLSLLGLAAFFSLSYVFIASIYLL  
 GEKLNHWKWGDMRQPRKKRK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_032744

**ORF Size:** 690 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_032744.4</a></u>
<b>RefSeq Size:</b>	1602 bp
<b>RefSeq ORF:</b>	693 bp
<b>Locus ID:</b>	84830
<b>UniProt ID:</b>	<u><a href="#">Q961Z2</a></u>
<b>Cytogenetics:</b>	6p24.1
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	Regulates the expression and the cell-associated anticoagulant activity of the inhibitor TFPI in endothelial cells (in vitro).[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG202724