

Product datasheet for **RG202792**

TBC1D7 (NM_016495) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: TBC1D7 (NM_016495) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: TBC1D7
Synonyms: MGCPH; PIG51; TBC7
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG202792 representing NM_016495
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACTGAGGACTCTCAGAGAACTTTCGTTTCAGTATATTATGAGAAAGTGGGGTTTCGTGGAGTTGAAG
 AAAAGAAATCATTAGAAATTCCTCTAAAAGATGACCGTCTGGATACTGAGAACTTTGTACTTTTAGTCA
 GAGGTTCCCTCTCCCGTCCATGTACCGTGCATTGGTATGGAAGGTGCTTCTAGGAATCTTGCCCTCCACAC
 CACGAGTCCCATGCCAAGGTGATGATGATCGTAAGGAGCAGTACTTGGATGTCCTTCATGCCCTGAAAG
 TCGTTCGCTTTGTTAGTATGCCACACCTCAGGCTGAAGTCTATCTCCGATGTATCAGTGGAGTCTGG
 GAAGTTACCTCGAAGTCCCTCTTTTCCACTGGAGCCAGATGATGAAGTGTTCCTTGCCATAGCTAAAGCC
 ATGGAGGAAATGGTGAAGATAGTGTGCGACTGTTACTGGATCACCCGACGCTTTGTGAACCAATTAATA
 CCAAGTACCGGGATTCCTTGCCCCAGTTGCCAAAAGCGTTTGAACAATACTTGAATCTGGAAGATGGCAG
 ACTGCTGACTCATCTGAGGATGTGTTCCGCGGCCGCCCCAACTTCTTATGATCTCTGGTTCAAGAGGTGC
 TTTGCGGGATGTTGCCTGAATCCAGTTTACAGAGGGTTGGGATAAAGTTGTGAGTGGATCCTGTAAAGA
 TCCTAGTTTTGTAGCTGCGAAATTTTAACTTTAAAATAAAGTTATGGCACTGAACAGTGCAGA
 GAAGATAACAAAGTTTCTGAAAAATATCCCCAGGACAGCTCAGACGCGATCGTGAGCAAGGCCATTGAC
 TTGTGGCACAACACTGTGGGACCCCGTCCATTCAAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG202792 representing NM_016495
Red=Cloning site Green=Tags(s)

MTEDSQRNFRSVYYEKVGFGRGVEEKKSLEILLKDDRLDTEKLCTFSQRFPLPSMYRALVWKVLLGILPPH
 HESHAKVMMYRKEQYLDVLHALKVVRFVSDATPQAEVYL RMYQLESGKLRSPSPFLEPDDEVFLAIKA
 MEEMVEDSVDCYWITRRFVNQLNTKYRDSLPLPKAFEQYLNLEDGRLLTHLRMCSAAPKLPYDLWFKRC
 FAGCLPESSLQRVWDKVVSGSCKILVFVAVEILLTFKIKVMALNSAEKITKFLENIPQDSSDAIVSKAID
 LWHKHCCTPVHSS

TRTRPLE - GFP Tag - V

Restriction Sites:


SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_016495

ORF Size: 879 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016495.6
RefSeq Size:	1322 bp
RefSeq ORF:	882 bp
Locus ID:	51256
UniProt ID:	Q9P0N9
Cytogenetics:	6p24.1
Gene Summary:	This gene encodes a member of the TBC-domain containing protein family. The encoded protein functions as a subunit of the tuberous sclerosis TSC1-TSC2 complex which plays a role in the regulation of cellular growth and differentiation. Mutations in this gene have been associated with autosomal recessive megalencephaly. Alternative splicing results in multiple transcript variants. Naturally occurring readthrough transcription occurs between this locus and downstream LOC100130357. [provided by RefSeq, Jan 2016]