

Product datasheet for **RG202820**

C11orf2 (VPS51) (NM_013265) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C11orf2 (VPS51) (NM_013265) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C11orf2
Synonyms:	ANG2; ANG3; C11orf2; C11orf3; FFR; PCH13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RG202820 representing NM_013265
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGGCAGCTGCCCGGGCCTAGCCCGGGTCTGGACCTGGGGACTCCCCAGAAGGGCCGAGG
 GGGAGGCTCCGGAGCGTCGGCGGAAGGCGCACGGGATGCTGAAGCTTTACTACGGCCTCTCGGAAGGGGA
 GGCGGGGACGCCCGCGGGGGCCGACCCCTGGACCCGACTGATCTGAACGGGGCGCACTTCGACCCG
 GAAGTTTACCTAGACAAGCTGCGTAGAGAGTGCCTCTGGCCAGTTGATGGACAGTGAGACGGACATGG
 TGCGGCAGATCCGGGCTCTAGACAGCGACATGCAGACCCTGGTCTATGAGAACTACAACAAGTTCATCTC
 AGCCACAGACACCATCCGGAAGATGAAGAACGATTTCCGGAAGATGGAGGATGAGATGGACCGGCTGGCC
 ACCAACATGGCAGTGATCACCGACTTCAGCGCTCGCATCAGCGCCACGCTGCAGGACCGCCACGAGCGCA
 TCACCAAGCTGGCAGGGTCCACGCGCTGCTGCGGAAGCTGCAGTTCCTCTTTGAGCTGCCCTCGCGCT
 CACCAAGTGCCTGGAAGTGGCGCCTATGGGCAGGCGGTGCGCTACCAGGGCCGCGCAGGCCGTGCTG
 CAGCAGTACCAACACCTGCCCTCGTCCGCGCCATCCAGGACGACTGCCAGGTATCACGGCCCCCTGG
 CCCAGCAGCTGCGGCAGCGCTTTAGGGAGGGCGGCTCAGGCGCCCCGGAGCAGGCAGAGTGCCTGGAGCT
 GCTGCTGGCCCTGGGCGAGCCTGCGGAGGAGCTGTGCGAGGAGTTCCTGGCGCACGCCCGCGCCGGCTG
 GAGAAGGAGCTGAGAACTGGAGGCCGAGCTGGGGCCCTCACCTCCGGCTCCCGACGTGTTAGAGTTCA
 CCGACCATGGAGGCAGTGGCTTCGTGGGCGGCCTCGCCAGGTGGCGGGCCCTACCAGGAGCTGTTTGC
 GGCCAGGGCCAGCAGGTGCCGAGAAGCTGGCGGCCTTCGCCGGCAGCTGGGCAGCCGCTATTTTGGC
 CTGGTGGAGCGGGCTGGCGCAGGAGCAGGGTGGTGGTGAACACTACTGCTGGTGCGGCGCTGGACC
 GCTTCCACCGCGCTTGGCGGCTCCCGGGCCCTGCTGGCCGCTGCCGGGCTCGCAGAGCTGCCACGGA
 GATCGTGAACGAGTGGCCCGCAGCGCCTGGGCCACCACCTGCAGGGTCTCCGGGCGGCCTTCCTGGGC
 TGCTGACAGACGTCGCGCAGGCGCTGGCAGCACCTCGCGTGGCTGGGAAGGAGGGCCCTGGCCTGGCCG
 AGTTGCTGGCAAATGTGGCCAGCTCCATCCTGAGCCACATTAAGGCCTCTCTGGCAGCAGTGACCTTTT
 CACCGCAAAGAGGTGTCCTTCTCAACAAGCCCTACTCCGGGGTGGTTCGACGTCAGGGTGTCCGT
 GAGGGCTCATCGTGGGCTTCGTCCACTCTATGTGCCAGACGGCTCAGAGCTTCTGCGACAGCCCTGGG
 AGAAGGGGGTGCCACACCACCTGCCCTGCTCCTGCTGCTCTCCCGCCTCTGCCTGGACTACGAGACGGC
 CACCATCTCTACATCCTCACTCTCACTGATGAACAGTTTCTGGTGCAGGATCAGTTCACAGTGACGCC
 GTGAGCACGCTGTGTGCAGAGGCCAGGGAACGGCGCGGGCTGCTGACCCACTACGTGAAGGTGCAGG
 GCCTGGTCATATCACAGATGCTGCGCAAGAGCGTGGAGACTCGCGACTGGCTCAGCACTCTGGAGCCCCG
 GAATGTGCGGGCCGTCATGAAGCGGGTGGTGGAGGATACCACCGCATCGACGTGCAGGTGGGGCTCCTG
 TACGAAGAGGGTGTTCGCAAGGCCCAGAGCAGCGACTCCAGCAAGAGGACTTTCTCCGTGTACAGCAGCT
 CTCGGCAGCAGGGCCGCTACGCCCCAGCTATACCCCAAGTGCCTCGATGGACACCAACCTTTGAGCAA
 TATCCAGAAGCTATTCTCTGAACGATTTGATGTGTTTCCAGCCCTGTGGAGTTCAACAAGGTGTCCGTGCTG
 ACCGGCATCATCAAGATCAGCCTGAAGACGCTGCTGGAGTGTGTGCGGCTGCGCACCTTTGGGCGCTTCG
 GGCTGCAGCAGGTGCAAGTGGACTGCCACTTTCTGCAGCTCTACCTGTGGCGTTTTGTGGCCGACGAAGA
 ACTCGTCACTTGCTGCTGGACGAAGTGGTGGCCTCTGCTGCCCTGCGCTGCCAGACCCCTGTGCCATG
 GAGCCAGTGTGGTTGAGGTCATCTGCGAGCGCGC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG202820 representing NM_013265
 Red=Cloning site Green=Tags(s)

MAAAAAAGSPGSGPGDSEPEGEAPERRRKAHGMLKLYYGLSEGEAAGRPAGPDPDPTDLNGAHFDP
 EVYLDKLRRECPALQMDSETDMVRQIRALDSDMQTLVYENYNKFI SATDTIRKMKNDFRKMEDMDRLA
 TNMAVITDF SARI SATLQDRHERITKLAGVHALLRKLQFLFELPSRLTKCVELGAYGQAVRYQGRAQAVL
 QQYQHLPSFRAIQDDCQVITARLAQQLRQRFREGGSGAPEQAECVELLLALGEPAEELCEEFLAHARGRL
 EKELRNLEAELGPSPAPDVLEF TDHGSGGFVGGLCQVAAAAYQELFAAQGPAGAEKLAAFARQLGSRFYFA
 LVERRLAQEQQGGDNSLLVRLDRFHRRLRAPGALLAAAGLADAATE IVERVARERLGHHLQGLRAAFLG
 CLTDVRQALAAPRVAGKEGPGLAELLANVASSILSHIKASLAAVHLFTAKEVSFSNKPYPFRGEFCSQGVV
 EGLIVGFVHSMCQTAQSFCDSPGEKGGATPALLLLL SRLCLDYETATISYILTLTDEQFLVQDQFPVTP
 VSTLCAEARETARRLL THYVKVQGLVISQMLRKS VETRDWLSTLEPRNVRVAVMKRVVETDTAIDVQVGLL
 YEEGVRKAQSSDSSKRTFSVYSSSRQQGRYAPSYTPSAPMDTNLLSNIQKLF SERIDVFSPEFVFNKVSVL
 TGIKISLKTLLCEVRLRTFGRFGLQQVQVDCHFLLQLYLWRFVADEELVHLLLDEVVASAALRCPDPVPM
 EPSVVEVICERG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_013265

ORF Size: 2346 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013265.4](#)

RefSeq Size: 2546 bp

RefSeq ORF: 2349 bp


Locus ID: 738

UniProt ID: [Q9UID3](#)

Cytogenetics: 11q13.1

Gene Summary: This gene encodes a member of the vacuolar protein sorting-associated protein 51 family. The encoded protein is a component of the Golgi-associated retrograde protein complex which acts as a tethering factor for carriers in retrograde transport from the early and late endosomes to the trans-Golgi network. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012]

Product images:



Circular map for RG202820