

Product datasheet for **RG202830**

MRPS30 (NM_016640) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPS30 (NM_016640) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MRPS30
Synonyms:	MRP-S30; PAP; PDCD9; S30mt
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG202830 representing NM_016640 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCGGCCAGGTGTTGGAGGCCTTTGCTACGCGTCCGAGGCTTTCATTGCACACCGCGGCTAATG
CCGCGGCCACGGCTACAGAAACGACCTGCCAAGACGTCGCGGCGACCCCGTCGCGCGGTACCGGCCGAT
TGTGGCCTCCATGACAGCCGACAGCAAGCTGCACGGCTGCGGCGGATCGAGCGCTGGCAGGCGACGGTG
CACGCTGCGGAGTCGGTAGACGAGAAGCTGCGAATCCTACCAAGATGCAGTTTATGAAGTACATGGTTT
ACCGCGAGACCTTCGCGCTGAATGCCGACCGCTGGTACCAGTACTTCACCAAGACCGTGTTCTGTCTGGG
TCTGCCGCCGCCCGGAGCGGAGCCGAGCCGAGCCGAAACCGAACCTGAACCTGCGCTGGACCTCGCG
GCGCTGCGTGCCTGCGCTGCGACTGCCTGCTGCAGGAGCACTTCTACCTGCGGCGCAGGCGGCGCTGC
ACCGTTACGAGGAGAGCGAGGTCAATCTTTGCCCTTCTCGATCAGCTGGTGTCAACCTCGTGGGCT
CCTCAGCCACACAACCGGCCCTGGCCGCTGCCGCCCTCGATTATAGATGCCAGTTTCATTTTACTGG
GTGCGTGGTGAAGAAATTATTCCTCGTGGTCATCGAAGAGGTGAATTGATGACTTGCATACCATAGATAG
ATGATAAACCAACAACCGAGATTCGAATATCCAAGCAACTCGCAGAGTTTGTGCCATTGGATTATCTGT
TCCTATAGAAATCCCCACTATAAAATGTAAACGAGACAACTTCCATTATCAAACGGCAGTATGAAAC
CACATATTTGTTGGCTCAAAACTGCAGATCCTTGCTGTTACGGTCACACCCAGTTTCATCTGTTACCTG
ACAAATTAAGAAGGGAAGGCTTTTGAGACAAACTGTGCTGATCAGATAGAAAGTTGTTTTAGAGCTAA
TGCTATTGCAAGCCTTTTGTCTGGACTGGAGCAAGCTATGTATCAAGGATTCTGGAGTGAAGCAGAT
GTTACTCGACCTTTTGTCTCCAGGCTGTGATCAGATGGAAAATACTTTCTTTTCTGCTACCAGC
TAAATACTTTGGCACTGACTACACAAGCTGATCAAAATAACCCTCGTAAAAATATATGTTGGGTACACA
AAGTAAGCCTCTTATGAAACAATTGAGGATAATGATGTGAAAGGTTTAAATGATGATGTTCTACTTCAG
ATAGTTCACTTTCTACTGAATAGACCAAAAGAAGAAAATCACAGCTGTTGAAAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG202830 representing NM_016640
 Red=Cloning site Green=Tags(s)

MAAARCWRPLLRGPRLSLHTAANAATATETTCQDVAATPVARYPPIVASMTADSKAARLRRIERWQATV
 HAAESVDEKLRIITKMQFMKYMVYPQTFALNADRWYQYFTKTVFLSGLPPPPAEPEPEPEPEPALDLA
 ALRAVACDCLLQEHFYLRRRRRVHRYEESEVISLPFLDQLVSTLVGLLSPHPALAAAALDYRCPVHFYW
 VRGEEIIPRGHRRGRIDDLRYQIDDKPNNQIRISKQLAEFVPLDYSVPPIEPTIKCKPKLPLFKRQYEN
 HIFVGSKTADPCCYGHTQFHLLPDKLRRERLLRQNCADQIEVVFRANAIASLFAWTGAQAMYQGFWSEAD
 VTRPFVSQAVITDGKYFSFFCYQLNTLALTTQADQNNPRKNICWGTQSKPLYETIEDNDVKGFNDVLLQ
 IVHFLNRPKEEKSQLEN

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_016640

ORF Size: 1317 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016640.4](#)

RefSeq Size: 1686 bp

RefSeq ORF: 1320 bp


Locus ID: 10884

UniProt ID: [Q9NP92](#)

Cytogenetics: 5p12

Gene Summary: Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that is similar to the chicken pro-apoptotic protein p52. Transcript variants using alternative promoters or polyA sites have been mentioned in the literature but the complete description of these sequences is not available. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG202830