

Product datasheet for **RG203159**

TRIP15 (COPS2) (NM_004236) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIP15 (COPS2) (NM_004236) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	COPS2
Synonyms:	ALIEN; CSN2; SGN2; TRIP15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG203159 representing NM_004236 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGACATGGAGGATGATTTTCATGTGCGATGATGAGGAGGACTACGACCTGGAATACTCTGAAGATA
GTAACCTCCGAGCCAAATGTGGATTTGGAAAATCAGTACTATAATTCCAAAGCATTAAAAGAAGATGACCC
AAAAGCGCATTAAAGCAGTTTCCAAAAGGTTTTGGAACCTGAAGGTGAAAAGGAGAATGGGGATTTAAA
GCACTGAAACAAATGATTAAGATTAACCTCAAGTTGACAACTTCCAGAAATGATGAATAGATATAAGC
AGCTATTGACCTATATTCGGAGTGCAGTCACAAGAAATTATTCTGAAAAATCCATTAATTCTATTCTTGA
TTATATCTCTACTTCTAAACAGATGGATTTACTGCAGGAATTCTATGAAACAACACTGGAAGCTTTGAAA
GATGCTAAGAATGATAGACTGTGGTTTAAAGACAAACAAAGCTTGAAAATTATATTTAGAACGAGAGG
AATATGAAAAGCTTCAAAAAATTTACGCCAGTTACATCAGTCGTGCCAGACTGATGATGGAGAAGATGA
TCTGAAAAAGGTACACAGTTATTAGAAATATATGCTTTGGAAATTCAAATGTACACAGCACAGAAAAAT
AACAAAAACTTAAAGCACTCTATGAACAGTCACCTCACATCAAGTCTGCCATCCCTCATCCACTGATTA
TGGGAGTTATCAGAGAATGTGGTGGTAAAATGCACTTGAGGGAAGGTGAATTTGAAAAGGCACACACTGA
TTTTTTTGAAGCCTTCAAGAATTATGATGAATCTGGAAGTCCAAGACGAACCACTTGCTTAAAATATTTG
GTCTTAGCAAATATGCTTATGAAATCGGGAATAAATCCATTTGACTCACAGGAGGCAAGCCGTACAAAA
ATGATCCAGAAATTTAGCAATGACGAATTTAGTAAGTGCCTATCAGAATAATGACATCACTGAATTTGA
AAAGATTTCAAAAACAAATCACAGCAACATCATGGATGATCCTTTATAAGAGAACACATTGAAGAGCTT
TTGCGAAACATCAGAACACAAGTCTTATAAAAATTAATTAAGCCTTACACAAGAATACATATTCCTTTTA
TTTCTAAGGAGTTAAACATAGATGTAGCTGATGTGGAGAGCTTGGTGCAGTGCATATTGGATAACAC
TATTCATGGCCGAATTGATCAAGTCAACCAACTCCTTGAAGTGGATCATCAGAAGAGGGGTGGTGCACGA
TATACTGACTAGATAAATGGACCAACCAACTAAATCTCTCAACCAGGCTGTAGTCAGTAACTGGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG203159 representing NM_004236
 Red=Cloning site Green=Tags(s)


MSDMEDDFMCDDEEDYDLEYSSEDSNSEPNVDLENQYYNSKALKEDDPKAALSSFQKVLLEGEKGEWGFK
 ALKQMIKINFKL TNFPMMNRYKQLLTYIRSAVTRNYSEKSINSILDYISTSKQMDLLQEFYETTLEALK
 DAKNDRLWFKTNTKLGKLYLEREEYGLKQILRQLHQSCQTDGDDLLKGTQLLEIYALEIQMYTAQKN
 NKKLKALYEQSLHIKSAIPHPLIMGVIRECGGKMHLEGEFEKAHTDFFFAFKNYDESGSPRRTTCLKYL
 VLNMLMKSGINPFDSQEAKPKYKNDPEILAMTNLVSAYQNDITEFEKILKTNHNSIMDDPFIREHIEEL
 LRNIRTQVLIKLIKPYTRIHIPFISKELNIDVADVESLLVQCILDNTIHGRIDQVNQLLELDHQKRGGAR
 YTALDKWTNQLNSLNQAVVSKLA

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_004236

ORF Size:	1329 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_004236.4
RefSeq Size:	1951 bp
RefSeq ORF:	1332 bp
Locus ID:	9318
UniProt ID:	P61201
Cytogenetics:	15q21.1
Domains:	PCI
Protein Families:	Druggable Genome, Transcription Factors
Gene Summary:	Essential component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. The complex is also involved in phosphorylation of p53/TP53, c-jun/JUN, IκappaBα/NFKBIA, ITPK1 and IRF8/ICSBP, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively. Involved in early stage of neuronal differentiation via its interaction with NIF3L1.[UniProtKB/Swiss-Prot Function]