

Product datasheet for **RG203280**

UTP23 (NM_032334) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	UTP23
Synonyms:	C8orf53
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG203280 representing NM_032334
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGAAGATCACAAGGCAGAAACATGCCAAGAAGCATCTTGGCTTCTCCGCAACAACCTCGGAGTCCGCG
AGCCGTACCAGATCCTGCTGGACGGCACCTTCTGTGAGCGCGCGCTGCGGGCCGCATCCAGCTGCGGGA
GCAGCTGCCCGCTACCTCATGGGGGAGACGCAGCTGTGCACCACAAGATGTGTGTTAAAGAGCTAGAA
ACATTGGGAAAGGACTTATATGGGGCAAACTGATTGCACAAAAATGCCAAGTTCGAAATTGCTCATT
TCAAGAATGCAGTGAGTGGATCAGAATGTCTGCTTCCATGGTTGAAGAGGGAATCCTCATCATTATTT
TGTGGCAACACAGGATCAGAATTTGTCTGTGAAAGTAAAAAGAAGCCTGGAGTTCCTCTCATGTTTATT
ATTGAGAACTATGGTTTGGACAAACCTTCTCCAAAACAATTGCCTTTGTAAGCAGTGGAGTCAG
GTCAGCTTGTCTCAGTGCATGAGAAAGAAAGTATCAACATCTCAAGAGGAACAGGGTTTGTGAAAAA
CACTGAACAGAGTAGAAGAAAAAGCGCAAGAAATAAGTGGTCCCAATCCTCTTAGTTGTTTGAAGAAA
AAGAAAAAGGCACCGGACACACAATCATCTGCTTCTGAAAAGAGAAAAAGAAAAAGAATTCCGAACA
GATCTAACCCAAAAGTACTTTCTGAGAAGCAGAATGCAGAAGGAGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG203280 representing NM_032334
Red=Cloning site Green=Tags(s)

MKITRQKHAKKHLGFFRNNFGVREPYQILLDGTFCQAALRGRIQLREQLPRYLMGETQLCTTRCVLKELE
TLGKDLYGAKLIAQKCQVRNCPHFKNVAVSGSECLLSMVEEGNPHHYFVATQDQNLVKVKKKPGVPLMFI
IQNTMVLDPKSPKTI AFVKAVESGQLVSVHEKESIKHLKEEQGLVKNTEQSRRKKRKKISGPNPLSCLKK
KKKAPDTQSSASEKRRKRKRIRNRSNPKVLSEKQNAEGE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_032334

ORF Size: 747 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_032334.1](#), [NP_115710.1](#)

RefSeq Size: 998 bp

RefSeq ORF: 750 bp

Locus ID: 84294


UniProt ID: [Q9BRU9](#)

Cytogenetics: 8q24.11

Domains: DUF652

Gene Summary: Involved in rRNA-processing and ribosome biogenesis. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG203280