

Product datasheet for **RG203317**

ATP6J (ATP6V1G1) (NM_004888) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ATP6J (ATP6V1G1) (NM_004888) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: ATP6J
Synonyms: ATP6G; ATP6G1; ATP6GL; ATP6J; Vma10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG203317 representing NM_004888
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTAGTCAGTCTCAGGGGATTGAGCAGCTGCTGCAGGCCGAGAAGCGGGCAGCCGAGAAGGTGTCCG
 AGGCCCGCAAAGAAAGAACCGGAGGCTGAAGCAGGCCAAAGAAGAAGCTCAGGCTGAAATTGAACAGTA
 CCGCCTGCAGAGGGAGAAAGAATTCAAGGCCAAGGAAGCTGCGGCATTGGGATCCCGTGGCAGTTGCAGC
 ACTGAAGTGGAGAAGGAGACCCAGGAGAAGATGACCATCCTCCAGACATACTCCGGCAGAACAGGGATG
 AAGTCTTGGACAACCTCTTGGCTTTTGTCTGTGACATTCGCCAGAAATCCATGAAAACCTACCGCATAAA
 TGGA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG203317 representing NM_004888
 Red=Cloning site Green=Tags(s)
 MASQSQGIQQLLQAEKRAAEKVSEARKRKNRRLKQAKEEAQAEIEQYRLQREKEFKAKEAAALGSRGSCS
 TEVEKETQEKMTILQTYFRQNRDEVLDNLLAFVCDIRPEIHENYRING

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_004888

ORF Size: 354 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_004888.4](#)

RefSeq Size: 1110 bp

RefSeq ORF: 357 bp

Locus ID: 9550
UniProt ID: [O75348](#)
Cytogenetics: 9q32
Domains: V-ATPase_G
Protein Pathways: Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection
Gene Summary: This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A, three B, and two G subunits, as well as a C, D, E, F, and H subunit. The V1 domain contains the ATP catalytic site. The protein encoded by this gene is one of three V1 domain G subunit proteins. Pseudogenes of this gene have been characterized. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG203317