

Product datasheet for **RG203355**

FLJ23584 (NM_024588) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: FLJ23584
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG203355 representing NM_024588
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAATCCTGTCCCCACTGGGGTGAAGTCTTCTGCTGGTGGTGGGAGGGAGAGCATCTGGCCAGCC
AGGGCACAACCCCTGCCAGGGATCATAGAGTGGGGATCAGTCCAGCCTCCAGCAGGCCCAACCTGAATC
CTGGAGAAGACGACAAAGGGACAAAGGGGTGGATCCAGAGAAGGCCCCAGCCTGACTCGGCAGTCCCAA
AACCCCTCCGTCTCTGACAGCTCCCTTGGGCATGCCCTCTGCCTGCTCCTGTCTCCCGTGTGGCCAGCCC
CGGAAGCTGCCATCATTCTAGCGGGTCTCCACCGCCCTCACTGTCTGCCCAAGGGGACAGGCCTCAA
GAAGAGCAAACGACTGCTCCTGGAGTCCCTCATGCGGAGGAGGATTGCACACCTGAAGTGGGGTCTTCCC
CGGCGGATCCTGGAGTCTATTTCTGTCTTAAGTCTTAGGATCTTGCTCATTGACCCCTGCTGGGGCGA
GGCTCTCTGGACTGAACACAGGCCAGGAGCTCCAAGCCAGCAGGAAAGGTATTGTGAGGCCCAAGGCTC
CCCACCGGCCTTAAGTCCCAGAGAGGTTCCAGAGGGTTCAGCGCCAGACAGAAAAGCTCGAACTT
CCTATACAAGCCAGAGCTCTGGAGAGGAACAGACCGCACATGTCAGAGCCCATTAAGCATTTCATCCAG
CC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG203355 representing NM_024588
Red=Cloning site Green=Tags(s)


MNPVPHWGEVFLLVGGEHGLASQGTTPARDHRVGISPASQQAQPESWRRRQRDKGVDPEKAPSLTRQSQ
NPPSLTAPLGMPACSLPCGPAPEAAIILAGPPTALTVLPKGTGLKKSKRLLESIMRRRIAHKWLPL
RRILESYFLNFLGSCSLTAGARLSGLNTGQELQAQQERYCEAQSPPLKSPERFQVRQRPDRKSSKL
PIQARALERNRPHMSEPIKHFHPA

TRTRPLE - GFP Tag - V



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_024588

ORF Size: 702 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_024588.3, NP_078864.1
RefSeq Size: 4992 bp
RefSeq ORF: 704 bp
Locus ID: 79640
Cytogenetics: 22q13.2

Product images:


Circular map for RG203355