

Product datasheet for **RG203404**

ZNF 556 (ZNF556) (NM_024967) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	ZNF 556
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG203404 representing NM_024967
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGACACAGTGGTCTTTGAAGACGTGGTTGTGGATTTACGCTGGAGGAGTGGGCCTTGCTGAATCCTG
CTCAGAGAAAACCTACAGAGATGTCATGCTGGAGACCTTCAAGCACCTGGCCTCAGTAGATAATGAGGC
TCAGCTTAAAGCCAGTGGGTCTATTTCTCAGCAGGATACCTTCTGGAGAAAAATTATCCCTCAAACAGAAA
ATAGAAAAGTTACAAGAAAGAATATATGGGCCTCCCTTTTAGGAAAAAATTGGGAAGAACATAGCGTTA
AAGACAAGCACAAACCAAGGAGAGACATTTGAGAAATCCAAGGGTGGAGAGACCATGTAAAAGCAGTAA
AGGTAATAAACGTGGAAGAACCTTCAGAAAGACTCGAAATTGTAATCGTCATCTGCGCAAGAATTGTTGT
ACTAGTGTAAGACGGTACGAATGCAGTCAGTGTGGAAAACTTTCACCCATTCCTCATCCCTGATAAGGC
ACAAAAGAGCTCACTCTGGACAAAAATTATATAATGTAAGGAATGTGGGAAAGCCTTCAGTCGCCCTTC
CTACCTACAGACGCATGAGAAAACCTCAGTGGAGAGAAACCCTATGCCTGTCAATCTTGCGGGAAGACA
TTTCTTCGTTCCCACTCTCTCACTGAACATGTAAGGACTCACACTGGAGAGAAACCCTACGAATGTGGGC
AGTGTGGGAAAGGCTTCAGTTGTCCCAATCCTTTCGCGCACATGTGATGATGCACGCCGAGGAGAGACC
GTATGAGTGCAAGCACTGTGGGAAAGCCTTCAGGTGTCAGAAATCCTTTCGAGTCCATATGATCATGCAC
GCCGGAGGGAGACCGTATGAGTGCAAGCAGTGTGGGAAAGCCTACTGCTGGGCAACATCCTTCAACGAC
ACGTGAGAATTACAACGGGGAGAAACCCTATAAGTGTGGAAAAATCGGGAAAGCATTCGGTTGGCCCTC
ATCCTTACACAAACACGCGAGAACGATGCTAAAAAGAAACCTGTGAGTGGGGGACGCGTGGGAAAGTCT
TCCGCGAGGCCTCGCCCTCCACAGATGTCAAATCACAACTAGAGAGAAAGTCTATAAATGTGAAACGT
GTGGGAAAACGTATGGTTGGTCTCATCTTACACAAACATGAGAGAAAGCACACTGGGGAGAAACCTGT
AAATGCAGCCAGTGTGGGAAAACCTTCAGGCGGGCTTTGCTCTTCAAAAATGTAAGAAGCAGATTGGA
CAGAAAGCCAGTAAATGCGAAAAATGTGGGAAAGCTTTCAGTTGTCCCAAGCCTTTCAGGTCATGTGA
GAAGTCACACAGGAAAGAAATCCTGTACATCTAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG203404 representing NM_024967

Red=Cloning site Green=Tags(s)


MDTVVFEDVVVDFTLEEWALLNPAQRKLYRDVMLETFKHLASVDNEAQLKASGSISQQDTSGEKLSLKQK
IEKFTRKNIWASLLGKNWEEHSVKDKHNTKERHLRNPRVERPCKSSKGNKRGRTFRKTRNCNRHLRKNCC
TSVRRYECSCQCKLFTSSSLIRHKRAHSGQKLYKCKECGKAFSRPSYLQTHEKTHSGEKPYACQSCGKT
FLRSHSLTEHVRTHTEGKPYECGQCGKGFSCPKSFAHVMMHAGGRPYECKHCGKAFRCQKSFVHMIMH
AGGRPYECKQCGKAYCWATSFQRHVRIHNGEKPYKCGKCGKAFGWPSLHKHARTHAKKPPVSGGSVGKS
SARPRPSTDVKSQTREKVYKCETCGKTYGWSSSLHKHERKHTGEKPVNAASVGKPSGGLCSSKNVRTQIG
QKPSKCEKCGKAFSCPKAFQGHVRSHTGKKSCTSK

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_024967

ORF Size:

1365 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_024967.1](#), [NP_079243.1](#)

RefSeq Size: 1632 bp

RefSeq ORF: 1371 bp

Locus ID: 80032


UniProt ID: [Q9HAH1](#)

Cytogenetics: 19p13.3

Protein Families: Transcription Factors

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG203404