

Product datasheet for **RG203419**

Destrin (DSTN) (NM_006870) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Destrin (DSTN) (NM_006870) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: Destrin
Synonyms: ACTDP; ADF; bA462D18.2; HEL32
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG203419 representing NM_006870
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTCAGGAGTGCAAGTAGCTGATGAAGTATGTCGCATTTTTATGACATGAAAGTTCGTAATGCT
 CCACACCAGAAGAAATCAAGAAAAGAAAGAAGGCTGTCATTTTTGTCTCAGTGCAGACAAAAAGTGCAT
 CATTGTAGAAGAAGGCAAAGAGATCTTGGTTGGAGATGTTGGTGAACCATAACTGATCCTTCAAGCAT
 TTTGTGGGAATGCTTCTGAAAAAGATTGTCGCTATGCTTTGTATGATGCAAGCTTTGAAACAAAAGAAAT
 CCAGAAAAGAAGAGTTGATGTTTTTTTTGTGGCACCAGAACTAGCACCTCTGAAAAGTAAAATGATCTA
 TGCAAGCTCCAAGGATGCAATTAAGAAAGAAATTTCAAGGCATAAACATGAATGTCAAGCAAATGGACCA
 GAAGATCTCAATCGGGCTTGTATTGCTGAAAAGTTAGGTGGATCCTTAATTGTAGCCTTTGAAGGATGCC
 CTGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


Protein Sequence: >RG203419 representing NM_006870
 Red=Cloning site Green=Tags(s)

MASGVQVADEVCRIFYDMKVRKCSTPEEIKKRKKAVIFCLSDKCKIIVEEGKEILVGDVGVITIDPFKH
 FVGMLEPKDCRYALYDASFETKESRKEELMFFLWAPELAPLKSMMIYASSKDAIKKKFQGIKHECQANGP
 EDLNRACIAEKLGGSLIVAFEGCPV

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_006870

ORF Size: 495 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006870.4](#)

RefSeq Size: 1614 bp


RefSeq ORF: 498 bp

Locus ID: 11034

UniProt ID: [P60981](#)

Cytogenetics:	20p12.1
Domains:	ADF
Gene Summary:	The product of this gene belongs to the actin-binding proteins ADF family. This family of proteins is responsible for enhancing the turnover rate of actin in vivo. This gene encodes the actin depolymerizing protein that severs actin filaments (F-actin) and binds to actin monomers (G-actin). Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG203419