

Product datasheet for **RG203424**

MRPS22 (NM_020191) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPS22 (NM_020191) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	MRPS22
Synonyms:	C3orf5; COXPD5; GIBT; GK002; MRP-S22; ODG7; RPMS22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG203424 representing NM_020191 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGCCCTCGGAACAACTGTATTGCTGTGGAGCCTCTTGAGGAGTTCTCCGGCGTGGAACGGGTCT
 GTTCCGGGCTCGAATCCAGCCCTGGCACGGTGGCCTGCTCAACCGCTACCTTGCTCTTCGAGATGGG
 GCTGCCACGCCCGGTTTCAGCTCCGAGGCCGAGAATCTGGTAGCCAGAGACCAAGAAGCCTACATTT
 ATGGATGAGGAAGTTCAAAGCACTACGAAAAATGACAGGCTTGAAGTGCAGAAGACTTTTAAGCCAG
 CTATACAAGAACTGAAGCCACCAACCTATAAGCTAATGACTCAGGCACAGTTGGAAGAGGCTACAAGACA
 GGCAGTTGAGGCAGCTAAAGTACGATTAATAATGCCACCAGTTCTGGAAGAGCGAGTACCAATAAATGAT
 GTGTTAGCTGAAGATAAGATTTTGAAGGAACAGAAACAACCAAAATATGTGTTTACTGATATATCATATA
 GCATACCACACCGGGAGCGTTTTATTGTCGTGAGAGAACCAAGTGGCACACTACGCAAAGCCTCTTGGGA
 AGAACGGGACCGAATGATACAAGTTTATTTCCAAAAGAAGGTCGTAAAAATTTGACACCAATAATTTTC
 AAGGAAGAAAATCTTAGGACTATGTATAGCCAGGACAGGCATGTTGATGTCCTCAATCTCTGCTTTGCC
 AGTTTGAGCCAGATTCACAGAGTATATCAAGGTTTCATCACAAGACCTATGAAGATATAGATAAACGTGG
 AAAATATGACCTTTTACGTTCAACAAGATACTTTGGTGAATGGTGTGGTATTTTGTAAATAATAAAAG
 ATTGATGGTTTGCTGATTGACCAGATTCAGAGAGATTTAATCGATGATGCAACCAACTTGGTCCAGCTGT
 ATCACGTGCTCCATCCAGATGGCCAGTCGGCTCAAGGGGCAAGGATCAGGCTGCTGAGGGAATAAATTT
 AATCAAGGTCTTTGCAAAAACAGAAGCACAGAAGGGAGCCTATATAGAACTAACACTGCAGACTTATCAA
 GAAGCACTCAGTCGCCATTCTGCAGCTTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >RG203424 representing NM_020191
 Red=Cloning site Green=Tags(s)


MAPLGTTVLLWSLLRSPGVERVCFRARIQPWHGGLLQPLPCSFEMGLPRRRFSSEAAESGSPETKKPTF
 MDEEVQSILTKMTGLNLQKTFKPAIQELKPPTYKLMTQAQLEEATRQAVEAAKVRMKMPVLEERVPIIND
 VLAEDKILEGTETTKYVFTDISYSIPHRERFIVREPSGTLRKASWEERDRMIQVYFPKEGRKILTPIIIF
 KEENLRMTYSQDRHVDVLNLCFAQFEPDSTEYIKVHHKTYEDIDKRGKYDLLRSTRYFGGMWYFVNNKK
 IDGLLIDQIQRDLDIDATNLVQLYHVLHPDGQSAQGAQKQAAEGINLIKVFATKTEAQKGAYIELTLQTYQ
 EALSRSAAAS

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_020191

ORF Size: 1080 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_020191.2](#), [NP_064576.1](#)

RefSeq Size: 1155 bp

RefSeq ORF: 1083 bp


Locus ID: 56945

UniProt ID: [P82650](#)

Cytogenetics: 3q23

Gene Summary: Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that does not seem to have a counterpart in prokaryotic and fungal-mitochondrial ribosomes. This gene lies telomeric of and is transcribed in the opposite direction from the forkhead box L2 gene. A pseudogene corresponding to this gene is found on chromosome Xq. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RG203424