

Product datasheet for **RG203495**

FAM3A (NM_021806) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM3A (NM_021806) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FAM3A
Synonyms:	2.19; DLD; DXS560S; XAP-7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG203495 representing NM_021806 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGTTGGCAGGCCCTCTCCGCATTGTGGTCTAGTCGTCAGTGTGGGTGTCACATGGATCGTGGTCA
GCATCCTCCTGGGTGGGCTGGCAGTGGCTTTCCTCGCATCCAGCAACTTTCACCAGTCCAGAGAGCTC
GGTGACTGCAGCGCCACGGCCAGGAAGTACAAGTGTGGCCTGCCCCAGCCGTGCTCTGAGGAGCACCTG
GCCTTCCGCGTGGTCAGCGGGCCGCCAACGTATTGGGCCAAGATCTGCCTCGAGGACAAGATGCTGA
TGAGCAGCGTCAAGGACAACGTGGGCCGCGGGCTGAACATCGCCCTGGTGAACGGGTGTCAGCGCGAGCT
CATCGAGGCCCGGCCTTTGACATGTGGGCCGGAGATGTCAACGACCTGTTGAAGTTTATTCGGCCACTG
CACGAAGGCACCCTGGTGTTCGTGGCATCTACGACGACCCAGCCACCAAGATGAATGAAGAGACCAGAA
AGCTCTTCAGTGAGCTGGGCAGCAGGAACGCCAAGGAGCTGGCCTTCCGGGACAGCTGGGTGTTTGTTCGG
GGCCAAGGGTGTGCAGAACAGAGCCCTTTGAGCAGCACGTGAAGAACAGTAAGCACAGCAACAAGTAC
GAAGGCTGGCCGAGGCGCTGGAGATGGAAGGCTGTATCCCGGAGAAGCACGGCCAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG203495 representing NM_021806
Red=Cloning site Green=Tags(s)

MRLAGPLRIVVLVVSIVGVTWIVVSIILLGGPGSGFPRIQQLFTSPSSVTAAPRARKYKCGLPQPCPEEHL
 AFRVVSQAANVIGPKICLEDKMLMSSVKDNVGRGLNIALVNGVSGELIEARAFDMWAGDVNDLLKFI
 RPLHEGTLVFVASYDDPATKMNEETRKLSELGSRNAKELAFRDSWVFGAKGVQNKSPFEQHVKNKSHSNKY
 EGWPEALEMEGCIPRRSTAS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021806

ORF Size: 690 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021806.4](#)

RefSeq Size: 1606 bp

RefSeq ORF: 693 bp

Locus ID: 60343


UniProt ID: [P98173](#)

Cytogenetics: Xq28

Protein Families: Secreted Protein, Transmembrane

Gene Summary: This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]

Product images:



Circular map for RG203495