



AGTCCCTGTGTGACCTGCGCCTCTCCCTCATTTCCCCACACGGCGGCTCTTGACCAGGCAGTGGGGG  
 TCGGGTGACCAGTATGGGACCTGAGGTGGTGTGCAGGCTGTGCCTTTGAAATTGATGGCTCTGAGGAG  
 ACTCTGGATTTCCACGGAGCTGGCTGTGCCTGTATCCGAGACCATGTTAGGAAACGCGACTTGTT  
 TTTTACCACCTACTTCTTGGCCCTGGCTAACACCTGAAGAGCAAAGCCATGGACCTGGCTCAGGCAGG  
 CAGCACAGTGAATCTAAGATCTACGACACTCCAGTGGCAGATGTGGACACTCCTGCCTGGTCTCTGC  
 ACAAGGCTACAGATGTGGCCATCTCTTCAAAGGGCTGGCACGGACGCTGGGCATGGCCATCAGCGAGC  
 GTCCAGACCTGAGGGTACCGTGTGCCAGGCCCTGCGCACCTCATACCAAGGGCTGCCAGGCAGAGGC  
 TGACCGTGTGAAGTGAGTCGCTTTGCCAAGAACTTTCTGCCGATCCTCTCAACCTGTATGGCAGCCC  
 GTGGCAGCCGGGACACTCCAGCCCCTCGCCGGCTGTGCTGGAACCATCAGAACTTACCTCACCATCA  
 CTGACACTCAGTTGGTGAACAGTCTCCTGAAAAAGCCAGTGAGAAGGTGCTCGACCCTGCCAGCTCTGA  
 AAGCTATACCTCCACCATCCGGCCCTACCTAGAGAGCAAAGGCCACGGGGTGCAGAAAGAGGCTACCGAG  
 TGCTGGAGGAGGTGTGTCCAGTCTCAGGGCCCCGGGCCCTCTTCTGTCAGAGCCACCTGGAGGACCT  
 GAAGAAGACACTGCTGGACTCGCTGCGGAGCACCTCCTACCCGCCAAGAGGCCCTTTGAAGTGCCTC  
 CTACACATCGTGAGGAAGCTCTCAGCTGAACACAAGGAGTTCATCACTGCCCTCATCCAGAGGTGATCC  
 TGTGCACCAAGGAGGTGTGCGTGGGCGCACGGAAGACGCTTTTGCACTGCTCGTGGAGATGGGCCATGC  
 TTTCTAAGGTTTGGCTCGAACCAGGAAGAGGCCCTGCAGTGTACCTCGTCTGATCTACCTGGCCTG  
 GTGGGCGCGGTGACCATGGTCAGCTGCAGCATCCTGGCCCTGACCCACCTCCTTTTCGAGTTTAAAGTCT  
 TGATGGGGACCACTACAGTGGAGCAGCTGCTGGAGAATGTGTGCTGCTTCTGGCCTCCGACCCGTTGA  
 CGTGGTCAAGTCTGCACTGGGCTTCATCAAGGTGGCAGTGACTGTATGGACGTGGCGCACCTGGCCAAA  
 CATGTGCAGCTGGTGTGAAGCCATTGGGAAGCTTTTCAAGTACATGCGGCGCACTTCCGCATGAAGC  
 TTCGGAACCTGTTACCAAGTTCATCCGCAAGTTTGGATTTGAGCTGGTAAAAGGCTGTTGCCGAGGA  
 GTACCACAGAGTCTGGTCAACATCCGAAAGCTGAGGCCCGGGCCAAGAGGCACCGAGCCCTGAGCCAG  
 GCTGCCGTGGAGGAGGAAGAAGAGGAGGAGGAGGAGGAGGCCCGCCAGGGCAAAGGTGACAGCATTG  
 AGGAGATTTTAGCTGACTCAGAGGACGAGGAGGACAATGAGGAGGAGGAAGAAGCCGAGGCAAGGAGCA  
 GCGGAAGCTGGCAGCAGAGGAGCGGGCATGGCTGAAAGAGGGCGGTGGGGACGAGCCCTCACTTC  
 CTGGATCCCAAGGTGGCCCAACGAGTCTTGCCACGCAGCCAGGGCCAGGCCGGGGCAGGAAGAAGGACC  
 ACAGCTTCAAGGTGAGCGCCGATGGCCGGCTGATCATAAAGGAGGAGGCAGACGGCAACAAGATGGAGGA  
 AGAGGAAGGTGCCAAAGGCGAAGATGAAGAGATGGCTGACCAATGGAAGATGTGATCATCAGGAATAAA  
 AAGCACCAGAAGCTCAAGCACCAGAAAGAGGCTGAGGAGGAGGAGTGGAGATACCCCTCAGTACCAAG  
 CTGGAGGCTCTGGCATTATCGCCCTGTGGCCAAGAAGGCTATGCCTGGGGCTGAATACAAGCCAAGAA  
 AGCAAAAGGTGATGTGAAGAAGAAAGGCCGGCCGGATCCCTATGCCTACATCCCCCTCAACAGAAGCAAG  
 CTCAACCGCAGGAAGAAGATGAAGCTGCAGGGACAGTTCAAAGGCTGGTGAAGGCTGCCCGGCGAGGTT  
 CCCAGGTGGGACAAAAACCGCAGAAAGGATCGTCGACCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>RG203502 representing NM\_015179

Red=Cloning site Green=Tags(s)


MGRSGKLP SGVS AKLKR WKKGHSSDSNP AICRHRQAARS RFFSRPSGRSDLTVDAVKLHNELQSGSLRLG  
KSEAPETPMEEAEVLTEKSSGTFLSGLSDCTNVTF SKVQRFWESNSAAHKEICAVLA AVEVIRSQGG  
KETET EYFAALMTTMEAVESPE SLAAVAYLLNLVLKRVSPVL IKKFSDTSKAFMDIMSAQASSGSTSVL  
RWVLSCLATLLRKQDLEAWGYPVTLQVYHGLLSFTVHPKPKIRKAAQHGVCSVLKGSEFMFEKAPAHHPA  
AISTAKFCIQEIEKSGGSKEATTTLHMLTLLKDLLPCFPEGLVKSCSETLLRVMTLSHVLVTACAMQAFH  
SLFHARPGLSTLSAELNAQIITALYDYVPSENDLQPLLAWLKVMEKAHINLVRLQWDLGLGHLPRFFGTA  
VTCLLSPHSQVLTAATQSLKEILKECVAPHMADIGSVTSSASGPAQSVAKMFRAVEEGLTYKFHAAWSSV  
LQLLCVFFEACGRQAHPVMRKCLQSLCDLRLSPHPHTAALDQAVGA AVTSMGPEVVLQAVPLEIDGSEE  
TLDFFRSWLLPVIRDHVQETRLGFFTTYFLPLANTLKSAMD LAQAGSTVESKIYDTLQWQMTLLPGFC  
TRPTDVAISFKGLARTLGMAISERPDLRVTVQC ALRTLITKGCQAEADRAEVS RFAKNFLPILFNL YGQP  
VAAGDTPAPRRAVLETIRTYLTITDTQLVNSLLEKASEKVLDPASSDFTRL SVLDLVVALAPCADEAAIS  
KLYSTIRPYLESKAHGVQKKAYRVLEEVCA SPQPGALFVQSHLEDLKKTLDSL RSTSSPAKRPRLKCL  
LHIVRKL SAEHKEFITALIPEVILCTKEVSVGARKNAFALLVEMGHAF LRFGSNQEEALQCYLVL IYPGL  
VGAVTMVSCSILALTHLLFEFKGLMGTSTVEQLLE NVCLLLASRTDRVVK SALGFIK VAVTVM DVHLAK  
HVQLVME AIGKLSDDMRHRFRM KLRNLFTK F IRKFGE L VKRLLPEEYHRVLVNI RKA EARA KRHRALSQ  
AAVEEEEEEEEEEP AQKGDSIEE ILADSEDEEDNEEEERSGKEQRKLARQ RSRAWLKEGGGDEPLNF  
LDPKVAQRVLATQPGPGRGRKKDHSFKVSADGRLI IREEADGNKMEEEEGAKGEDEEMADPMEDVIIRNK  
KHQKLKHQKEAEEEELEIPPQYQAGGSGIHRPVAKKAMPGA EYKAKKAKG DVKKKGRDPYAYIPLNRSK  
LNRRKKMKLQGQFKGLVKAARRGSQVGHKNRRKDRRP

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:


NM\_015179

ORF Size:

3891 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>                                                                      |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <a href="#">NM_015179.2</a> , <a href="#">NP_055994.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq Size:</b>           | 4387 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 3894 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 23223                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <a href="#">Q5JTH9</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 10q24.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

Product images:



Circular map for RG203502