

Product datasheet for **RG203682**

C4orf29 (ABHD18) (NM_001039717) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	C4orf29
Synonyms:	C4orf29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG203682 representing NM_001039717
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGTGAGCAAGTTAGATATTCTATACCGGAGACTTCTCTAACAAAACCTTTTATCAGAGGATGGG
GAAGGCCAGAAGATCTCAAAAGACTCTTTGAATTCAGAAAGATGATTGGAATCGGGAAGATGCCAGAA
TCTGGTTTCAAGCGATTATCCAGTACACATTGATAAGATTGAAGAGCAATCAGATTGTAAGATCTTAGAT
GGACACTTTGTTTCCCCATGGCTCACTATGTGCCTGATATCATGCCAATTGAATCTGTTATTGCAAGGT
TCCAATTTATTGTGCCTAAAGAATGGAACAGCAAATATAGACCTGTATGCATTATCTTGCTGGAACAGG
AGATCATCATTACTGGAGGCGACGAACACTAATGGCCCGTCTATGATTAAAGAAGCCGAATGGCTTCT
TTATTGTTAGAAAACCTTATTATGGCTGCAGGAAACCAAGGACCAAGTAAGGTCCAGCTTAAAAAATG
TGTCCGACCTTTTGTGATGGGAGGAGCTCTTGTTTAGAATCTGCAGCTCTTGCAGTGGCTAGAGAG
GGAAGGTTACGGCCCTTTAGGAATGACTGGAATATCCATGGGAGGACACATGGCTTCCTTAGCGGTATCC
AACTGGCCTAAGCCCATGCCATTGATTCCATGCCTGTCTTGGTCCACAGCATCTGGGGTCTTCACTACGA
CAGATTCTTTCAAAATGGGACAAGAGTTTGTGAAACACTTCACTAGCAGTGCAGACAAGCTAACTAACCT
TAATCTGGTTTCCAGAACTTTAAATTTAGATATATCAAACCAAGTTGTATCCCCAAAACCTGCTGACTGC
CATAATTTAGCAAAACATCTGTCACTGCGACATCAGAAGGACTCTTATTGCAAGATACCTCTAAGATGA
AGCGCTTCAATCAACACTTTCAACCAACAAAAGTGGTTATACAAGTCGCAACCCTCAGTCATACCACCT
ACTTAGTAAGAACAAGCAGAAACAGTCTTCGAAAGAGTCTTTAATTTATGAAAGGAGTCATGGAT
GAATGTACTCATGTAGCAAATTTCTCAGTCCAGTTGATCCAAGCCTCATTATAGTGGTTCAAGCCAAAG
AAGATGCCTATATTCCACGAACAGGAGTTCGAAGTTTACAAGAAATTTGGCCTGGTTGTGAAATCCGATA
CTTAGAAGGGGTCAATTAGTGCTTATCTTTTAAACAAGGACTCTTCAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence:

>RG203682 representing NM_001039717

Red=Cloning site Green=Tags(s)


MGVSKLDILYRRLLLTKLFIKRWGRPEDLKRLFQFRKMIGNRERCQNLVSSDYPVHIDKIEEQSDCKILD
GHFVSPMAHYVPDIMPISVIAFQFIVPKEWNSKYRPVCIHLAGTGDHYYWRRRTLMARPMIKEARMAS
LLENPPYYGCRPKDQVRSSLKNVSDLFVMGGALVLESAALLHWLEREGYPLGMTGISMGHMASLAVS
NWPKPMPILIPCLSWSTASGVFTTTDSFKMGQEFVKHFTSSADKLTNLNLVSRTLNLDISNQVVSQKPADC
HNSSKTSVSATSEGLLLQDTSKMKRFNQTLSTNKSQYTSRNPQSYHLLSKEQSRNSLRKESLIFMKGVMD
ECTHVANFVSPVPSLIIVVQAKEDAYIPRTGVRSLQEIWPGCEIRYLEGGHISAYLFKQGLFR

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001039717

ORF Size:

1242 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001039717.2


RefSeq Size: 2200 bp

RefSeq ORF: 1245 bp

Locus ID: 80167

Cytogenetics: 4q28.2

Product images:



Circular map for RG203682