

Product datasheet for **RG203709**

EWSR1 (NM_005243) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EWSR1 (NM_005243) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	EWSR1
Synonyms:	bK984G1.4; EWS; EWS-FLI1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG203709 representing NM_005243
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGGCGTCCACGGATTACAGTACCTATAGCCAAGCTGCAGCGCAGCAGGGCTACAGTGCTTACACCGCCC
AGCCCACTCAAGGATATGCACAGACCACCCAGGCATATGGGCAACAAAGCTATGGAACCTATGGACAGCC
CACTGATGTCAGCTATACCCAGGCTCAGACCCTGCAACCTATGGGCAGACCCGCTATGCAACTTCTTAT
GGACAGCCTCCCCTGTTATACTACTCCAACCTGCCCCCAGGCATACAGCCAGCCTGTCCAGGGGTATG
GCACTGGTGCTTATGATACCACCCTGCTACAGTACCACCACCCAGGCCTCCTATGCAGCTCAGTCTGC
ATATGGCACTCAGCCTGCTTATCCAGCCTATGGGCAGCAGCCAGCAGCCACTGCACCTACAAGACCCGAG
GATGGAACAAGCCCACTGAGACTAGTCAACCTCAATCTAGCACAGGGGGTTACAACCAACCCAGCCTAG
GATATGGACAGAGTAACACAGTTATCCCAGGTACCTGGGAGCTACCCATGCAGCCAGTCACTGCACC
TCCATCCTACCCTCCTACCAGCTATTCTCTACACAGCCGACTAGTTATGATCAGAGCAGTTACTCTCAG
CAGAACACCTATGGGCAACCCAGCAGCTATGGACAGCAGAGTAGCTATGGTCAACAAAGCAGCTATGGGC
AGCAGCCTCCCCTAGTTACCCACCCAACTGGATCCTACAGCCAAGCTCCAAGTCAATATAGCCAACA
GAGCAGCAGCTACGGGCAGCAGAGTTTATTCCGACAGGACCACCCAGTAGCATGGGTGTTTATGGGCAG
GAGTCTGGAGGATTTCCGGACCAGGAGAGAACCCGAGCATGAGTGGCCCTGATAACCCGGGCAGGGGAA
GAGGGGGATTTGATCGTGGAGGCATGAGCAGAGGTGGCGGGGAGGAGGACGCGGTGGAATGGGCCTGG
AGAGCGAGGTGGCTTCAATAAGCCTGGTGGACCATGGATGAAGGACCAGATCTTGATCTAGGCCACCT
GTAGATCCAGATGAAGACTCTGACAACAGTGAATTTATGTACAAGGATTAATGACAGTGTGACTCTAG
ATGATCTGGCAGACTTCTTTAAGCAGTGTGGGGTTGTTAAGATGAACAAGAGAAGTGGCAACCCATGAT
CCACATCTACCTGGACAAGGAAACAGGAAAGCCAAAGGCGATGCCACAGTGTCCCTATGAAGACCCACCC
ACTGCCAAGGCTGCCGTGGAATGGTTTGTGGGAAAGATTTTCAAGGGAGCAAACCTTAAAGTCTCCCTTG
CTCGGAAGAAGCCTCCAATGAACAGTATGCGGGGTGGTCTGCCACCCCGTGGGGCAGAGGCATGCCACC
ACCACTCCGTGGAGGTCCAGGAGGCCAGGAGGTCTGGGGGACCCATGGGTGCGATGGGAGGCCGTGGA
GGAGATAGAGGAGGCTTCCCTCCAAGAGGACCCCGGGTTCGGAGGGAACCCCTCTGGAGGAGGAAACG
TCCAGCACCGAGCTGGAGACTGGCAGTGTCCCAATCCGGGTGTGGAAACCAGAACTTCGCTGGAGAAC
AGAGTGAACCAAGTGAAGGCCCAAGCCTGAAGGCTTCTCCCGCCACCCCTTCCGCCCCGGGTGGT
GATCGTGGCAGAGGTGGCCCTGGTGGCATGCGGGGAGGAAGAGGTGGCCATGGATCGTGGTGGTCCCG
GTGGAATGTTAGAGGTGGCCGTGGTGGAGACAGAGGTGGCTTCCGTGGTGGCCGGGCATGGACCAGG
TGCTTTGGTGGAGGAAGACGAGGTGGCCCTGGGGGGCCCTGGACCTTGTGGAACAGATGGGAGGA
AGAAGAGGAGGACGTGGAGGACCTGAAAAATGGATAAAGGCGAGCACCGTCAAGGAGCGCAGAGATCGGC
CCTAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG203709 representing NM_005243
 Red=Cloning site Green=Tags(s)

```

MASTDYSTYSQAAAQQGYSAITAQPTQGYAQTQAYGQQSYGTYGQPTDVSYTQAQTTATYGQTAYATSY
GQPPTGYTTPTAPQAYSQPVQGYGTGAYDTTATVTTTQASYAAQSAYGTQPAYPAYGQQPAATAPTRPQ
DGNKPTETSQPQSSTGGYNQPSLGYGQSNYSYPQVPGSYPMQPVTAAPPSPPTSYSTQPTSYDQSSYSQ
QNTYGPSSYGGQSSYGGQSSYGGQPPPTSYPPTSGYSYQAPSQYSQQSSSYGQQSSFRQDHPSSMGVYVQ
ESGGFSGPGENRSMSPDNRGRGRGGFDRGGMSRGGRRGGRRGGMGAGERGGFNKPGGPMDEGPLDLGPP
VDPDESDNSAIYVQGLNDSVTLDDLADFFKQCGVVKMNKRTGQPMIHIYLDKETGKPKGDATVSYEDPP
TAKAAVEWFDGKDFQGSKLKVS LARKKPPMNSMRGGLPPREGRGMPPLRGGPGGPGGGPMGRMGGRG
GDRGGFPFRGSRGNPSGGGNVQHRAGDWQCPNPGCGNQNF AWRTECNQCKAPKPEGFLPPPFPFPGG
DRGRGGPGGMRGGRLMDRGGPGGMFRGGRRGGDRGGFRGGRRGMDRGGFGGRRGGPGGPPGLMEQMGG
RRGRRGGPGKMDKGEHRQERRDRPY
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_005243

ORF Size: 1965 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005243.2 , NP_005234.1
RefSeq Size:	2644 bp
RefSeq ORF:	1971 bp
Locus ID:	2130
UniProt ID:	Q01844
Cytogenetics:	22q12.2
Domains:	RRM, zf-RanBP
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transcription Factors
Gene Summary:	<p>This gene encodes a multifunctional protein that is involved in various cellular processes, including gene expression, cell signaling, and RNA processing and transport. The protein includes an N-terminal transcriptional activation domain and a C-terminal RNA-binding domain. Chromosomal translocations between this gene and various genes encoding transcription factors result in the production of chimeric proteins that are involved in tumorigenesis. These chimeric proteins usually consist of the N-terminal transcriptional activation domain of this protein fused to the C-terminal DNA-binding domain of the transcription factor protein. Mutations in this gene, specifically a t(11;22)(q24;q12) translocation, are known to cause Ewing sarcoma as well as neuroectodermal and various other tumors. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 1 and 14. [provided by RefSeq, Jul 2009]</p>

Product images:



Circular map for RG203709