

Product datasheet for **RG203887**

C11orf54 (NM_014039) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C11orf54 (NM_014039) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C11orf54
Synonyms:	PTD012; PTOD012
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG203887 representing NM_014039 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGGCTTGTGCTGAGTTTTCTTTTCATGTACCAAGTCTTGAAGAGCTTGCTGGAGTTATGCAGAAGGGGT
TAAAAGATAACTTTGCTGATGTCCAGGTCTCTGTAGTTGATTGCCCTGATTTGACTAAGGAACCTTTAC
CTTTCCTGTAAGGCATCTGTGGGAAACTAGAATTGCAGAAGTTGGAGGTGTGCCTTACTTATTGCCT
CTTGTAACCAAAAAAAGTTTATGATCTGAATAAAATTGCAAAAGAAATCAAGCTGCCTGGAGCCTTTA
TTCTTGGAGCAGGAGCAGGTCCATTTAGACTCTCGGGTTCAATTCTGAGTTTATGCCAGTTATTCAGAC
AGAAAGTGAACACAAGCCTCCTGTAATGGAAGTTACTTTGCCCATGTGAACCCTGCAGATGGAGGGTGC
CTACTGGAGAAATACAGTGAGAAATGTCATGATTTTCAGTGTGCATTACTGGCTAATCTTTTTGCCAGTG
AAGGCCAACCTGGCAAGCCTGCAGAATTTCTTCTGCCCTTGAAGTCTGATGAAGAAGTGAATAAATG
GTTGCATTTTTATGAAATGAAAGCTCCTTTGGTTGTCTACCAGTTTTGTCTCCAGAGACCCAGGGTTT
GATTTGCGACTGGAGCACACTCATTTTTTTAGTCGTCATGGAGAAGGTGGACACTACCATTATGACACTA
CTCCAGATATAGTGAATATCTGGATACTTCTTACCTGCAGAGTTTCTCTATCGCATTGATCAACCAA
AGAGACGCATTCAATTGGCGGAGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_014039.4](#)

RefSeq Size: 2403 bp

RefSeq ORF: 798 bp


Locus ID: 28970

UniProt ID: [Q9H0W9](#)

Cytogenetics: 11q21

Gene Summary: Exhibits ester hydrolase activity on the substrate p-nitrophenyl acetate.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG203887