

Product datasheet for **RG203940**

GSG1L (NM_144675) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: GSG1L (NM_144675) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: GSG1L
Synonyms: PRO19651
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG203940 representing NM_144675
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGTCTCGAGCTCTTCCACTCCAGCAATGTCATCGACGGGCTCAAGCTCAATGCCTTCGCGGCTGTCT
TCACGGTGCTCTCAGGCCTCCTGGGAATGGTCGCCACATGATGTACACGCAGGTGTTCCAGGTACCCGT
GAGCCTCGGTCTGAGGACTGGAGACCCATTCTGGGACTACGGGTGGTCTTCTGCCTGGCGTGGGGC
TCCTTTACCTGCTGCATGGCAGCCTCTGTCACCACGCTCAACTCCTACACCAAGACGGTCATTGAGTTCC
GGCACAAGCGCAAGGTCTTTGAGCAGGGCTACCGGAAGAGCCGACCTTCATAGACCCTGAGGCCATCAA
GTACTTCCGGGAGAGGATGGAGAAGAGGGACGGGAGCGAGGAGGACTTCACTTAGACTGCCGCCACGAG
AGATACCTGCCCGACACACAGCCACACATGGCGGATTCCTGGCCCCGGAGCTCCGCACAGGAAGCACCAG
AGCTGAACCGACAGTGTGGGTCTTGGGCACTGGGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG203940 representing NM_144675
Red=Cloning site Green=Tags(s)

MCELFHSSNVIDGLKLNFAAVFTVLSGLLGMVAHMMYTQVFQVTVSLGPEDWRPHSWDYGWSFLAWG
SFTCCMAASVTTLNSYTKTVIEFRHKRVFEQGYREPTFIDPEAIKYFRERMEKRDGSEEDFHLDRCRHE
RYPARHQPHMADSWPRSSAQEAPELNRQCWVLGHWV

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_144675

ORF Size: 528 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144675.2](#)

RefSeq Size: 1188 bp

RefSeq ORF: 531 bp

Locus ID: 146395


UniProt ID: [Q6UXU4](#)

Cytogenetics: 16p12.1

Protein Families: Transmembrane

Gene Summary: As a component of the inner core of AMPAR complex, modifies AMPA receptor (AMPA) gating.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG203940