

Product datasheet for **RG204070**

SFT2D1 (NM_145169) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SFT2D1 (NM_145169) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: SFT2D1
Synonyms: C6orf83; pRGR1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG204070 representing NM_145169
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGAAGCTGCGGCGAGTCTGAGCGGCCAGGACGACGAGGAGCAGGGCCTGACTGCGCAGGTCCTGG
 ATGCCTCATCCCTTAGTTTCAACACCAGATTGAAATGGTTTGCCATCTGCTTCGTATGTGGCGTTTTCTT
 TTCTATTCTTGGAACTGGATTGCTGTGGCTTCCGGCGGCATAAAGCTTTTTGCAGTGTTTTATACCCTC
 GGCAATCTTGCTGCGTTAGCCAGTACATGCTTTTTAATGGACCTGTGAAGCAACTGAAGAAAATGTTTG
 AAGCAACAAGATTGCTTGAACAATTGTTATGCTTTTGTGTTTCATATTTACCCTGTGTGCTGCTCTTTG
 GTGGCATAAAGAAGGACTGGCTGTGTTATTCTGCATATTGCAGTTCTTGTCAATGACCTGGTATAGCCTG
 TCATACATCCCATATGCAAGGGATGCAGTTATTAATGCTGTTCTTCTCTCTAAGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG204070 representing NM_145169
 Red=Cloning site Green=Tags(s)

MEKLRRVLSGQDDEEQGLTAQVLDASSLSFNTRLKWFALCFVCGVFFSILGTGLLWLPGGIKLFAVFYTL
 GNLAALASTCFLMGPVKQLKMFATRLLATIVMLLCFIFTLCAALWWHKGLAVLFCILQFLSMTWYSL
 SYIPYARDAVIKCCSSLLS

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_145169

ORF Size: 477 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_145169.1, NP_660152.1</u>
RefSeq Size:	717 bp
RefSeq ORF:	480 bp
Locus ID:	113402
UniProt ID:	<u>Q8WV19</u>
Cytogenetics:	6q27
Protein Families:	Transmembrane
Gene Summary:	May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.[UniProtKB/Swiss-Prot Function]