

Product datasheet for **RG204366**

GLI4 (NM_138465) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	GLI4
Synonyms:	HKR4; ZNF928
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG204366 representing NM_138465
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCAGCCCTAGGGGACATTCAGGAGTCCCCTTCTGTCCCGTCCCCTGTCACTCTCTCATACCCGGGA
CACCTGGAACCCAGCACACGAGCCTCAGCTTCACTCCATGGGCATCAACATGGCTCCCCTGGCTCCAG
CCCTAAGGTGCTCTCCAGCCGTCCGACCTGGATCTCCAAGACGTAGAGGAAGTGGAGATCGGCAGAGAC
ACCTTCTGGCCGACTCCGAGCCCAAGCCGAGCAGGCTCCACGCTCTCCTGGCTCTCAGGCCCTGACG
AGGGGGCGGGCGGGCGCTGCGCAGCCTCCTGAGGAGCCTTCCCCGAGGGCCCGGTGCAGCGCCGGCTT
CGGGCCTGAATCCAGCGCGAGCGGCCGCGGGCCAGCCGCTGGGGCCGTCCCTTGCGCCAGCCGCGG
GGCGCCTGGCGGTGACGCTCGTGCAGCAAGCAGCGGCCGGGCCGAGGGTGCGCCGAGCGGGCTGCCG
AGCTGGGAGTCAACTTCGGTTCGGAGCCGGCAGGGCAGCGCGGGGGGCCAAGCCGACAGGTGCGAGGC
CTGCGCAAGAGTTTCAAGTATAACTCGCTGCTCCTGAAGCACCAGCGCATCCACACGGGCAGAAAGCCC
TACGCTGCCACGAGTGCAGCAAGCGCTTCCGCGGCTGGTTCGGGCTTCATCCAGCACCACCGCATCCACA
CGGGCGAGAAGCCCTACGAGTGCAGCCAGTGCAGCCGCGCTTCAGCCACAGCTCGCACTTCACGCAGCA
CCTGCGCATCCACAACGGCGAGAAGCCCTACAAGTGCGGCGAGTGCAGCCAGGCCTTCAGCCAGAGCTCC
AACCTGGTGCGCCACCAAGCGGCTGCACAGGGTGAGAAGCCCTACGCTGCAGCCAGTGCAGCAAGCCCT
TCATCTGGAGCTCCGTGCTCATCGAGCACCAGCGCATCCACACTGGCGAGAAGCCCTACGAGTGTCCGA
CTGCGGCAAGCCCTCCGCGGCGCTCGCACTTCTCCGGCACCTGCGGACCCACAGGGCGAGAAGCCC
TTCGCGTGTGGCGCTGCGGCAAGCCCTTCGGCCAGAGCTCCAGCTCATCCAGCACCAGCGGTGCACT
ACCGCGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAA



Protein Sequence: >RG204366 representing NM_138465
Red=Cloning site Green=Tags(s)

MAALGDIQESPSVPSVLSPPGTPGTQHHEPQLHLHGHQHGSPGSSPKVLSPQSDLDLQDVVEEVEIGRD
TFWPDSEPKPEQAPRSPGSQAPDEGAGGALRSLRLRRARCSAGFGPESSAERPAGQPPGAVPCAQPR
GAWRVTLVQQAAGPEGAPERAAELGVNFGSRQGSARGAKPHRCEACGKSFKYNSLLKHKQRIHTGEKP
YACHECGKRFRGWSGFIQHHRIHTGEKPYECGQCGRAFSSSHFTQHLRIHNGEKPYKCGECGQAFSSQSS
NLVRHQRLHTGEKPYACSGCGKAFIWSSVLIHQRIHTGEKPYECSDCGKAFGRSHFFRHLRTHHTGEKP
FACGACGKAFGQSSQLIQHQRVHYRE

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_138465


ORF Size: 1128 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_138465.4
RefSeq Size:	1363 bp
RefSeq ORF:	1131 bp
Locus ID:	2738
UniProt ID:	P10075
Cytogenetics:	8q24.3

Product images:



Circular map for RG204366