

Product datasheet for **RG204566**

FOXR2 (NM_198451) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOXR2 (NM_198451) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	FOXR2
Synonyms:	FOXN6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG204566 representing NM_198451 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGACTTAAACTAAAAGACTGTGAATTTTGGTATAGTCTCCATGGCCAGGTCCCAGGGCTGCTGGACT
GGGACATGAGGAATGAGTTATTTCTGCCTTGACACAGACCAGTGCTCTTAGCTGAGCAAATCCTTGC
CAAATACAGAGTCGGAGTAATGAAGCCCCAGAAATGCCTCAGAAGAGGAGACCCAGTCTGATGGAGAT
GGTCTCCCTGTGAACCAATCTGTGGATGTGGTGGACCCCAATATCCTGTGCCCCCTTGGCAGCCAGG
AGGCCCAAAGCCAGTGAAAAGAGGATCTGACAAACATTTCTCCTTTCCCTCAGCCCCACAAAAGA
CGAAGGGTCTAACTGCTCAGAGGACAAAGTGGTAGAGTCTCTGCCATCTTCTCCAGTGAGCAGTCTCCT
TTACAGAAGCAGGGTATCCATTCCCCAGTGACTTTGAGCTCACAGAAGAGGAGGCTGAGGAACACAGCG
ACAACTCCCTCCAGTCCCCTGAAATGAAATGTTACCAGAGCCAGAACTATGGCAAATCAACAACCAAGA
GAAGTCCTGGCAAAGGCCCTCTCAATTGTAGCCACCTATTGCCCTAGCATTAAAGAAACAACCCAC
TGTGGCCTCAGTGTGAGGAGATCTACAATTTACCCGACAGCATTTCCTTTTCTGGACAGCTCCGG
ATGGCTGGAAGAGCACCATTATTACAACCTCTGCTTCTGGACAGCTTTGAGAAGGTGCCAGACAGCCT
TAAGGATGAAGATAATGCAAGACCTCGCTCTTGCTTTGGAAGCTCACTAAGGAGGGCACCGCCGCTTT
TGGGAGGAGACTCGTGTCTTAGCCTTTGCTCAAAGGGAGAGAATCCAAGAGTGCATGAGTCAGCCAGAGT
GTTGACCTCTCTTTGATCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >RG204566 representing NM_198451
 Red=Cloning site Green=Tags(s)

MDLKLKDCFWYSLHGQVPGLLDWMRNEFLPCTTDQCSLAEQILAKYRVGMKPPKRRPSPDGD
 GPPCEPNLWMWVDPNLCPLGSQEAPKPSGKEDLTNISFPFPQPPQKDEGSNCSEDKVVESLPSSSSEQSP
 LQKQGIHSPSDFELTEEEAEEPDDNSLQSPKCYQSQKLWQINNQEKSQRPLNCSHLIALALRNPH
 CGLSVQEIYNFTRQHFFFWTAPDGWKSTIHYNLCFLDSFEKVPDSLKDEDNARPRSLWKLTKEGHRRF
 WEETRVLAFAQRERIQECMSQPELLTSLFDL

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_198451

ORF Size: 933 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_198451.4](#)

RefSeq Size: 1954 bp

RefSeq ORF: 936 bp


Locus ID: 139628

UniProt ID: [Q6PJQ5](#)

Cytogenetics: Xp11.21

Protein Families: Transcription Factors

Product images:



Circular map for RG204566