

Product datasheet for **RG204755**

FAM122C (NM_138819) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Tag: TurboGFP

Symbol: FAM122C

Mammalian Cell Neomycin

Selection:

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG204755 representing NM_138819
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCACAGGAGAAAATGAACTAGGTTTCAAGTCGCTGCCGAGTTCCACTACCGCAGACGGCAACATTC
TGAGAAGAGTCAACAGTGCCCTTTGATCAATGGACTTGGTTTAAATTCACAGGTGTTGCAAGCTGACAT
GTTAAGAATTAGGACAAACAGAACATTTAGGAATCGACGCTCTCTGTTGTTGCCACCTCCTCCCTTT
CATGGTTCCATCAGCCGCCTTCATCAAAATCAAACAGGAAGAAGCCATGGATTTAATAAATAGAGAAACAA
TGCTGAATGGAAGCTACAAAGTGAGATACAGATAAGTCACCTCTTGGGAAGAAGGCTTGAACTGAATGA
CAATGGCTTACAGAAATCATCCTCTCTAAAGTGCATTGATTAACTCCAGTATCCTCAATGGCTTCTCC
ATCAAGAAGACTGGGAAGTCTACTTCCAGCTACTTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTAA

Protein Sequence: >RG204755 representing NM_138819
Red=Cloning site Green=Tags(s)


MAQEKMKLGFKSLPSSTADGNILRRVNSAPLINGLFNSQVLQADMLRIRTNRTTFRNRRSLLPPPPF
HGSISRLHQIKQEEAMDLINEITMSEWKLQSEIQISHSWEEGLKLNENGLQKSSSLKCIDLTPVSSMASS
IKKTGKSTSSYF

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_138819

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_138819.4
RefSeq Size: 1370 bp
RefSeq ORF: 459 bp
Locus ID: 159091
UniProt ID: Q6P4D5

Cytogenetics: Xq26.3

Product images:



Circular map for RG204755