

Product datasheet for **RG204977**

serum amyloid A2 (SAA2) (NM_030754) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: serum amyloid A2 (SAA2) (NM_030754) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: serum amyloid A2
Synonyms: SAA; SAA1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG204977 representing NM_030754
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTTCTCACGGGCTGGTTTTCTGCTCCTTGGTCTGAGTGTGAGCAGCCGAAGCTTCTTTTCGT
TCCTTGGCGAGGCTTTTGGTGGGCTCGGGACATGTGGAGAGCCTACTCTGACATGAGAGAAGCCAATTA
CATCGGCTCAGACAAATACTTCCATGCTCGGGGAATATGATGCTGCCAAAAGGGACCTGGGGGTGCC
TGGGCCGAGAAGTGATCAGCAATGCCAGAGAGAATATCCAGAGACTCACAGGCCATGGTCCGGAGGACT
CGCTGGCCGATCAGGCTGCCAATAAATGGGCGAGGAGTGGCAGAGACCCCAATCACTTCCGACCTGCTGG
CCTGCCTGAGAAATAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG204977 representing NM_030754
Red=Cloning site Green=Tags(s)
MKLLTGLVFCSLVLSVSSRSFFSFLGEAFDGDARMWRAYSMDREANYIGSDKYFHARGNYDAAKRGPGGA
WAAEVISNARENIQRLTGHGAEDSLADQAANKWGRSGRDPNHFRPAGLPEKY


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_030754

ORF Size: 366 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_030754.2](#), [NP_110381.1](#)

RefSeq Size: 526 bp

RefSeq ORF: 369 bp

Locus ID: 6289
UniProt ID: [P02735](#)
Cytogenetics: 11p15.1
Protein Families: Druggable Genome
Gene Summary: This gene encodes a member of the serum amyloid A family of apolipoproteins. The encoded preproprotein is proteolytically processed to generate the mature protein. This protein is a major acute phase protein that is highly expressed in response to inflammation and tissue injury. This protein also plays an important role in HDL metabolism and cholesterol homeostasis. High levels of this protein are associated with chronic inflammatory diseases including atherosclerosis, rheumatoid arthritis, Alzheimer's disease and Crohn's disease. This protein may also be a potential biomarker for certain tumors. Finally, antimicrobial activity against *S. aureus* and *E. coli* resides in the N-terminal portion of the mature protein. [provided by RefSeq, Jul 2020]

Product images:



Circular map for RG204977