

## Product datasheet for **RG204995**

### ETV3 (NM\_005240) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ETV3 (NM\_005240) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ETV3  
**Synonyms:** METS; PE-1; PE1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG204995 representing NM\_005240  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAAGCCGGCTGTAGCATCGTGGAAAAGCCAGAAGGAGGTGGAGGTATCAGTTTCCTGACTGGCCCT  
ACAAAACAGAGTCATCCCAGGCTCCCGGCAGATCCAGCTGTGGCACTTCATCCTGGAGCTGCTGCAGAA  
GGAAGAGTTCCGCCATGTCATCGCCTGGCAGCAGGGAGAGTACGGGGAATTTGTCATCAAGGATCCAGAT  
GAGGTGGCCCGCCTCTGGGGCCGCAGGAAATGCAAACCACAGATGAATTATGACAAGCTGAGCCGGGCC  
TCAGATACTATTACAACAAGAGGATCCTTCATAAAACAAAAGGGAAAAGATTTACCTATAAATTTAACTT  
CAACAAGCTGGTATGCCAACTACCCATTCATCAACATTCGGTCAAGTGGTAAGATACAACTCTTTTG  
GTAGGGAAT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >RG204995 representing NM\_005240  
Red=Cloning site Green=Tags(s)

MKAGCSIVEKPEGGGYQFPDWAYKTESSPSRQIQLWHFILELLQKEEFRHVIWQQGEYGEFVIKDPD  
EVARLWGRRKCKPQMNYDKLSRALRYYYNKRILHKTKGKRFTYKFNFNKLVMPNYPFINIRSSGKIQTLL  
VGN

**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_005240

**ORF Size:** 429 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_005240.2](#)

**RefSeq Size:** 1443 bp

**RefSeq ORF:** 432 bp

**Locus ID:** 2117


**UniProt ID:** [P41162](#)

**Cytogenetics:** 1q23.1

**Protein Families:** Transcription Factors

**Gene Summary:** Transcriptional repressor that contribute to growth arrest during terminal macrophage differentiation by repressing target genes involved in Ras-dependent proliferation. Represses MMP1 promoter activity.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for RG204995