

Product datasheet for RG205010

SDHC (NM_003001) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: SDHC (NM_003001) Human Tagged ORF Clone

Tag: TurboGFP

Symbol: SDHC

Synonyms: CYB560; CYBL; PGL3; QPS1; SDH3

Mammalian Cell

Selection:

Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide >RG205010 representing NM_003001

Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGGCTGCGCTGTTGCTGAGACACGTTGGTCGTCATTGCCTCCGAGCCCACTTTAGCCCTCAGCTCTGTA
TCAGAAATGCTGTTCCTTTGGGAACCACGGCCAAAGAAGAGATGGAGCGGTTCTGGAATAAGAATATAGG
TTCAAACCGTCCTCTGTCTCCCCACATTACTATCTACAGTTGGTCTCTTCCCATGGCGATGTCCATCTGC
CACCGTGGCACTGGTATTGCTTTGAGTGCAGGGGTCTCTTTTTTGGCATGTCGGCCCTGTTACTCCCTG
GGAACTTTGAGTCTTATTTGGAACTTGTGAAGTCCCTGTTCTGGGGCCAGCACTGATCCACACAGCTAA
GTTTGCACTTGTCTTCCCTCTCATGTATCATACCTGGAATGGGATCCGACACTTGATGTGGGACCTAGGA
AAAGGCCTGAAGATTCCCCAGCTATACCAGTCTGGAGTTGTCCTGGTTCTTACTGTGTTGTCCTCTA

TGGGGCTGGCAGCCATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG205010 representing NM_003001

Red=Cloning site Green=Tags(s)

MAALLLRHVGRHCLRAHFSPQLCIRNAVPLGTTAKEEMERFWNKNIGSNRPLSPHITIYSWSLPMAMSIC HRGTGIALSAGVSLFGMSALLLPGNFESYLELVKSLCLGPALIHTAKFALVFPLMYHTWNGIRHLMWDLG

 ${\sf KGLKIPQLYQSGVVVLVLTVLSSMGLAAM}$

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-Mlul



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200


CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Cloning Scheme:





ACCN: NM_003001

ORF Size: 507 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 003001.5</u>

RefSeq Size: 1315 bp
RefSeq ORF: 510 bp
Locus ID: 6391
UniProt ID: Q99643



Cytogenetics: 1q23.3

Domains: Sdh_cyt

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Alzheimer's disease, Citrate cycle (TCA cycle), Huntington's disease, Metabolic pathways,


Oxidative phosphorylation, Parkinson's disease

Gene Summary: This gene encodes one of four nuclear-encoded subunits that comprise succinate

dehydrogenase, also known as mitochondrial complex II, a key enzyme complex of the tricarboxylic acid cycle and aerobic respiratory chains of mitochondria. The encoded protein is one of two integral membrane proteins that anchor other subunits of the complex, which form the catalytic core, to the inner mitochondrial membrane. There are several related pseudogenes for this gene on different chromosomes. Mutations in this gene have been associated with paragangliomas. Alternatively spliced transcript variants have been

described. [provided by RefSeq, May 2013]

Product images:



Circular map for RG205010