

## Product datasheet for **RG205100**

### hnRNP F (HNRNPF) (NM\_004966) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	hnRNP F (HNRNPF) (NM_004966) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	hnRNP F
Synonyms:	HNRPF; mcs94-1; OK/SW-cl.23
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG205100 representing NM_004966 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATGCTGGGCCCTGAGGGAGGTGAAGGCTTTGTGGTCAAGCTCCGTGGCCTGCCCTGGTCTGCTCTG  
TTGAGGACGTGCAGAACTTCTCTCTGACTGCACGATTTCATGATGGGGCCGAGGTGCCATTTTCATCTA  
CACTAGAGAGGGCAGGCAGAGTGGTGGGCTTTGTTGAAGTGGATCAGAAGATGATGAAAAATGGCC  
CTGAAAAAAGACAGGGAAAGCATGGGACACCGGTACATTGAGGTGTTCCAGGTCCACAGAACCGAGATGG  
ATTGGGTGTTGAAGCACAGTGGTCCCAACAGTGCCGACAGCGCCAACGATGGCTTCGTGCGGCTTCGAGG  
ACTCCCATTTGGATGCACAAAGGAAGAAATTTGTCAGTTCCTCTCAGGGTTGAAAATGTGCCAAACGGG  
ATCACATTGCCTGTGGACCCGAAGCAAGATTACAGGGGAAGCGTTCGTGCAGTTTGCCTCGCAGGAGT  
TAGCTGAGAAGGCTCTAGGGAAACACAAGGAGAGGATAGGGCACAGGTACATTGAGGTGTTAAGAGCAG  
CCAGGAGGAAGTTAGGTCATACTCAGATCCCCCTCTGAAGTTCATGTCCGTGCAGCGGCCAGGGCCCTAT  
GACCGGCCCGGGACTGCCAGGAGGTACATTGGCATCGTGAAGCAGGCAGGCCTGGAAGGATGAGGCCTG  
GTGCCTACAGCACAGGCTACGGGGGCTACGAGGAGTACAGTGGCCTCAGTGATGGCTACGGCTTACCAC  
CGACCTGTTCCGGAGAGACCTCAGTACTGTCTCTCCGGAATGTATGACCACAGATACGGCGACAGTGAG  
TTCACAGTGCAGAGCACACAGCCACTGTGCCATGAGGGGCTGCCGTACAAGCGACCCGAGAACG  
ACATTTACAACCTTCTCTCTCTCAACCCTGTGAGAGTCCATATTGAGATTGGCCAGATGGAAGAGT  
GACGGGTGAAGCAGATGTTGAGTTTGCTACTCATGAAGAAGTGTGGCAGCTATGTCCAAAGACAGGGCC  
AATATGCAGCACAGATATATAGAATCTTCTTGAATTCACAACAGGGGCCAGCAATGGGGCGTATAGCA  
GCCAGGTGATGCAAGGCATGGGGGTGTCTGCTGCCAGGCCACTTACAGTGGCCTGGAGAGCCAGTCAGT  
GAGTGGCTGTTACGGGGCCGCTACAGTGGGCAGAACAGCATGGGTGGCTATGAC

**ACGGTACGGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG205100 representing NM\_004966  
 Red=Cloning site Green=Tags(s)

```

MMLGPEGGEGFVVKLRGLPWSCSVEDVQNFLSDCTIHDGAAGVHFITYTREGRQSGEAFVELGSEDDVKMA
LKKDRESMGHRYIEVFRSHRTEMDWLKHSGPNSADSANDGFVRLRGLPFGCTKEEIVQFFSGLEIVPNG
ITLPVDPEGKITGEAFVQFASQELA EKALGKHKERIGHRYIEVFKSSQEEVRSYSDPPLKFMSVQRPQPY
DRPGTARRYIGIVKQAGLERMRPGAYSTGYGGYEEYSGLSDGYGFTTDLFGRDLSYCLSGMYDHRYGDSE
FTVQSTTGHCVHMRGLPYKATENDIYNFFSPLNPVRVHIEIGPDGRVTGEADVEFATHEEAVAAMSKDRA
NMQHRYIELFLNSTTGASNGAYSSQVMQGMVSAQAATYSGLESQSVSGCYGAGYSGQNSMGGYD
  
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_004966

**ORF Size:** 1245 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_004966.2](#), [NP\\_004957.1](#)

**RefSeq Size:** 2591 bp

**RefSeq ORF:** 1248 bp

**Locus ID:** 3185


**UniProt ID:** [P52597](#)

**Cytogenetics:** 10q11.21

**Domains:** RRM

**Gene Summary:** This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has three repeats of quasi-RRM domains that bind to RNAs which have guanosine-rich sequences. This protein is very similar to the family member hnRPH. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]

### Product images:



Circular map for RG205100