

Product datasheet for **RG205277**

POMT1 (NM_001077365) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	POMT1 (NM_001077365) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	POMT1
Synonyms:	LGMD2K; LGMDR11; MDDGA1; MDDGB1; MDDGC1; RT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG205277 representing NM_001077365
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTGGGGATTTTTGAAGCGCCCTGTAGTGGTGACGGCTGACATCAACTTGAGCCTTGTGGCCCTGACTG
GGATGGGGTTACTGAGCCGGCTGTGGCGACTCACCTACCAGCGGGCTGTGGTTTTTGACGAAGTATATTA
TGGGCAGTACATCTCTTTTTACATGAAACAAATCTTCTTCTGGATGACAGTGGGCCGCCATTTGGCCAC
ATGGTGTGTCCTTGGGAGGTTATTTAGGAGGATTCGATGGCAATTTTTGTGGAACAGAATTGGAGCAG
AATACAGTAGCAACGTGCCTGTGTGGTCCCTGCGCCTGCTGCCAGCACTCGCGGGGGCTTGTGGTCCC
CATGGCCTACCAGATAGTGTGGAGCTCCACTTTTCTCATTGTGCCGCCATGGGAGCTGCTCTGTTGATG
CTTATCGAGAATGCTCTCATCACTCAGTCAAGGCTAATGCTTTTGAATCAGTGTAATATTTTTCAATC
TATTGGCCGTGTGTCTACCTGAAGTTCTTCACTGCCAAGAGCACAGCCCTTTTCTCTGAGCTGGTG
GTTCTGGCTAACACTGACAGGGTCTGTTGCTTCTGTGCAAGTGGCATCAAGTACATGGGTGTTCACG
TACGTGCTCGTGTGGTGTGACGCTGTCCATGCCTGGCACCTGCTTGGAGACCAGACTTTGTCCAATG
TCTGTGTGTTCTGTCACTTGCTCGCCGAGCAGTGGCTTTGCTGGTCAATCCCGGTGCTCCTGTACTTACT
GTTCTTCTACGTCCACTTGATTCTAGTCTTCCGCTCTGGGCCCCACGACCAAATCATGTCCAGTGCCTTC
CAGGCCAGCTTAGAGGGAGGACTAGCTCGGATCACCCAGGGTACGCCACTGGAGGTGGCCTTTGGGTCCC
AGGTCACTCTGAGGAACGTCTTTGGAAACCTGTGCCCTGCTGGCTTCAATCCCACCAGGACACCTACCC
CATGATATATGAGAACGGCCGAGGCAGCTCCCACCAGCAACAGGTGACCTGTTACCCCTTCAAAGACGTC
AATAACTGGTGGTGTAAAGGATCCCAGGAGGCACAGCTGGTGGTGGCAGCCCTCCGAGACCTGTGA
GGCAGCGGGACATGGTGCAGTGGTCCAGGCATGACCACCCGCTCCCTGAACACGCATGATGTTGCAGC
CCCCCTGAGCCCCATTACAGGAGGTCTCCTGCTACATTGACTATAACATCTCCATGCCCGCCAGAAC
CTCTGGAGACTGGAAATTGTGAACAGAGGATCTGACACAGACGTCTGGAAGACCATCCTCTCAGAGGTCC
GCTTTGTGCAGTGAACACTTCCGCTGTCTAAAGCTGAGCGGGGCTCACCTCCCTGACTGGGGGTATCG
GCAACTGGAGATCGTCGGGGAGAAGCTGTCCCGGGGCTACCACGGGAGCACGGTGTGGAACGTGGAGGAG
CACCGATACGGCGCAGCCAGGAGCAGAGGGAGCGGGAACGGGAGCTGCACTCACCTGCGCAGGTGGACG
TCAGCAGGAACCTCAGCTTCATGGCGAGATTCTCGGAGCTGCAGTGGAGGATGCTGGCGCTGAGAAGTGA
TGACTCGGAACACAAGTACAGCTCCAGCCACTGGAGTGGTACCCTGGACACCAATATTGCCTACTGG
CTGACCCAGGACCAGCGCTCAGATCCACCTACTTGGAAACATAGTGATCTGGGTTTCGGGCAGCCTCG
CTCTGGCCATCTACGCCCTGCTGTCTTGTGGTACCTGCTCCGACGGCGAAGAAATGTCCATGACCTCCC
TCAGGATGCCTGGCTGCGCTGGGTGCTGGCTGGGGCGCTGTGTGCCGGTGGCTGGGCAGTGAACCTC
CCGTTCTTCTGATGGAGAAGACACTTCTCTACTACCTACCTACCCGCACTCACCTTCCAATCCTTC
TGCTCCCTGTGGTCTGCAGCACATCAGCGACCACCTGTGCAGGTCCCAGCTCCAGAGGAGCATCTTCAG
CGCCCTGGTGGTGGCCTGGTACTCCTCCGCGTCCACGTGTCCAACACGCTGCGCCCACTCACCTACGGG
GACAAGTCACTCTCGCCACATGAACTCAAGGCCCTTCGCTGGAAGACAGCTGGGACATCTTGATCCGAA
AACAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >RG205277 representing NM_001077365
 Red=Cloning site Green=Tags(s)


MWGFLKRPVVVTADINLSLVALTGMGLLSRLWRLTYQRAVVFDEVYVYQYISFYMKQIFFLDDSGPPFGH
 MYLALGGYLGFDGNFLWNRIGAEYSSNPVWSLRLPALAGALSVPMAVQIVLELHFHSHCAAMGAALLM
 LIENALITQSRMLLESVLIFFNLLAVLSYLKFFNCQEHSPFSLSWFWLTLTGACSCAVGIKYMVFVT
 YVLVLGVAAVHAWHLLGDQTL SNVCVFCHLLARAVALLVIPVVL YLLFFVYVHLILVFRSGPHDQIMSSAF
 QASLEGGLARITQGQPLEVAFGSQVTLRNVFGKPVPCWLHSHQDTYPMIYENGRGSSHQQVTCYPFKDV
 NNWVVVKDPRRHQLVVSPPRPVRHGMVQLVHGTTTSLNTHDVAAPLSPHSQEVSCYIDYINISMPAQN
 LWRLEIVNRGSDTDVWTKILSEVRFVHVNTSAVLKLSGAHLPDWGYRQLEIVGEKLSRGYHGSTVWVVEE
 HRYGASQEQRERERELHSPAQVDVSRNLSFMARFSELQWRMLALRSDSEHKYSSSPLEWVTLDTNIAYW
 LHPRTSAQIHLLGNIVIWVSGSLALAIYALLSLWYLLRRRRNVHDL PQDAWL RWVLGALCAGGWAVNYL
 PFFLMEKTLFLYHYLPALTFQILLLPVVLQHISDHLCRSQRSLFISALVVAWYSSACHVNTLRPLTYG
 DKSLSPHELKALRWKDSWDILIRKH


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_001077365

ORF Size: 2175 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001077365.1](#), [NP_001070833.1](#)

RefSeq Size: 3080 bp

RefSeq ORF: 2178 bp

Locus ID: 10585

UniProt ID: [Q9Y6A1](#)

Cytogenetics: 9q34.13

Protein Families: Transmembrane

Protein Pathways: O-Mannosyl glycan biosynthesis

Gene Summary: The protein encoded by this gene is an O-mannosyltransferase that requires interaction with the product of the POMT2 gene for enzymatic function. The encoded protein is found in the membrane of the endoplasmic reticulum. Defects in this gene are a cause of Walker-Warburg syndrome (WWS) and limb-girdle muscular dystrophy type 2K (LGMD2K). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2008]