

Product datasheet for **RG205408**

CCER1 (NM_152638) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	CCER1
Synonyms:	C12orf12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG205408 representing NM_152638
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACTCAGACCCTCGACACAAGGGAAGACCCTCTGAACCTGGGCGCGGTGGCGGCGCGGCTGTGGCT
GTGGCTGGGCACACTCGGCCTCCTTGAGCTCCTGGTCGTCCTGCCATCGAAGGCGCCCGGTGCTCCAGC
GTACAATAGACCGCACCAGATATAGCCCCAAGACCGAGTATGGGCCCAAGGAAGCAGCCGAAGCAACAG
CACGGCCCGGGCTTTTGGTTCCAACACCCGTGTGCTCTAACTGGGGTGCTGGGAGGGCCCTGGCGCC
CACCCTCCAGGATTCTGGAATTCCTGCGCGGTGCAAGTGTTCGGGTGTATGGCTGCACCTCT
CTGCTTTGCTGCTCCTGCTGGAGCGGTCTGGAACCTGGCTGGTGAAGCCCCAGGCAGGAAG
AAGCGCTGGGCGCGCAGAGGCCGCGGCTGCGCCACCACCCTCGCCACTCCTACCGCGGAGCCCGCCAG
CGGATGTGAGCACGCTGCCGCGCCGTCAAGCTGTATGAGTGGAGAGAGCCTGGCATGCGAGCGCCGCC
CAACACCACCAATTATCATGAACAGATCTACGAGGACATGAGGCAGCAGGAGAAGGTGGAGCGTCAG
CAGGAGGCGCTGCGGGCGCAGAAGGCCACGGTGAGCGGCGAGGCTCCCGAGCCAGATCCTCCGAAACG
ACGCGCCCCCTGGCGGCAGCAAGGAAACCTGGGACTGCAGGAACTCTGTATGGCTTTGTGAGAATCC
CTCTTAGCATTAGTCCCAACCCAGAGGAAAACAGTCTCTTCCCCGCTGCTGGTGAAGAAGAGGAG
GAGAAGAAAAATGATGATGAGGAGGAGTATGACCAGGAGGTGTGTATGCAAAGGAGGCGAGCGAGGAGG
AAGAAGAGGTCGAAGATGAGGAGGAAGAGGTGCAAGATGAGGAGGAAGAAGAGGTGCAAGAGGCTGAATA
TGTGGAGGAGGAGAGGAGGAGCTGGAAGAGGAGGAGCTGGAAGAGGAAGAGGAGTCTGGAGGAGAAC
GAGCAGAGAGGGAAGAATTTCACTTGCCTCTGGAATGCCTTTATCAATCTTCGTAGAGGTGAAGAAA
AGAGAGAGAACTTTATAAGCTGCACCTTTTAAACCCAGAGCAGATAATCCCAAAGTGCCACAGGAATC
CTGTTCATGGCACAGGACTTTAACTGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG205408 representing NM_152638
Red=Cloning site Green=Tags(s)

MTQTLDTREDPLNLGGGGGGCGCGWAHSASLSSWSSCHRRRPGAPAYNRPHRYSPKTEYGPPRKQPKQQ
HGPWFQPPVCSNWGCWGGPWRPPPGFWKFCPVQVFRVYGLHPLCFCCSCWGSWNPGWVKPPGRK
KRWGRRGRGLRHHPRHSYPRSPADVSTLPRPVKLIEWREPGMRAPPNTTQFIMNQIYEDMRQKEKVERQ
QEALRAQKATVSGEASPARSSGNDAPPGSKETWGLQETLYGFVQNP SLAFSPNPEENQSLAPLLVEEEE
EKKNDDDEEYDQEVCDAKEASEEEEEVEDEEEVEDEEEVEEAEYVEEGEELEEELEEEVELEEN
EQRGEFHLPLEPLSIFVEAEKRENFISCTFLNPEQIIPKVPQESLFMAQDFNC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_152638

ORF Size: 1218 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_152638.4

RefSeq Size: 2977 bp


RefSeq ORF: 1221 bp

Locus ID: 196477

UniProt ID: Q8TC90

Cytogenetics: 12q21.33

Product images:



Circular map for RG205408