

Product datasheet for **RG205410**

CAPSL (NM_001042625) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CAPSL (NM_001042625) Human Tagged ORF Clone
Tag: TurboGFP
Symbol: CAPSL
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >RG205410 representing NM_001042625
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGATCCAGGCCAAGAAAAAGCTCACCACGGCCACCAACCCATTGAAAGACTCCGACTGCAGTGCC
TGGCCAGGGGCTCTGCTGGGATCAAAGGACTTGGCAGAGTGTTTGAATTATGGATGACGATAATAATCG
AACCCCTTGATTTTAAAGAATTTATGAAAGGGTTAAATGATTATGCTGTGGTCATGGAAAAAGAAGAGGTG
GAAGAAGCTTTCCGGAGGTTTGATAAAGATGGAATGGAACAATAGACTTCAATGAATTTCTTCACAT
TAAGACCTCCAATGTCCAGAGCCAGAAAAGAGGTAATCATGCAAGCTTTTAGAAAGTTAGACAAGACTGG
AGATGGTGTATATAACAATCGAAGACCTTCGTGAAGTATATAATGCAAAACACCAACCAAGTACCAGAAT
GGGGAATGGAGTGAGGAACAAGTATTTAGGAAATTTCTGGATAACTTTGATTACCCCTATGACAAAGATG
GATTGGTGACCCCTGAGGAGTTCATGAAGTACTATGCAGGTGTGAGCGCATCCATTGACACTGATGTGTA
CTTCATCATCATGATGAGAACCCTGGAAGCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG205410 representing NM_001042625
Red=Cloning site Green=Tags(s)

MAIQAKKKLTTATNPIERLRLQCLARGSAGIKGLGRVFRIMDDNNRRLDFKEFMKGLNDYAVVMEKEEV
EELFRFDKDGNGTIDFNEFLLTLRPPMSRARKEVIMQAFRKLDKTGDGVITIEDLREVYNKHPKYQN
GEWSEEQVFRKFLDNFSPYDKDGLVTPEEFMNYAGVSASIDTDVYFIIMMRTAWKL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001042625

ORF Size: 624 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042625.1](#), [NP_001036090.1](#)

RefSeq Size: 987 bp


RefSeq ORF: 627 bp

Locus ID: 133690

UniProt ID: [Q8WWF8](#)

Cytogenetics: 5p13.2

Product images:



Circular map for RG205410