

Product datasheet for **RG205478**

LYZL4 (NM_144634) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Tag: TurboGFP
Symbol: LYZL4
Synonyms: LYC4; LYZA
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG205478 representing NM_144634
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAAGGCATCCGTGGTTCTCTCCCTCCTTGGCTACCTGGTGGTCCAAGTGGTGCTTACATCTTGGGGC
GTTGCACAGTGGCTAAGAACTCCACGATGGAGGCCTGGATTATTTGAGGGCTATAGCCTTGAGAACTG
GGTGTGCCTGGCCTACTTCGAGAGCAAGTTCAACCCCATGGCCATCTACGAGAACACACGTGAGGGCTAC
ACTGGCTTTGGCCTCTTTCAGATGCGTGGCAGTGACTGGTGTGGCGACCATGGCAGGAACCGCTGCCATA
TGTCATGTTCCGCTTTACTGAATCCTAATTTAGAGAAGACAATTAATGTGCCAAGACCATTGTAAAAGG
AAAAGAAGGGATGGGAGCATGGCCACCTGGTCCCGTACTGCCAGTACTCCGATACCCTGGCACGGTGG
CTGGATGGTTGCAAGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG205478 representing NM_144634
Red=Cloning site Green=Tags(s)


MKASVVLSSLGYLVVPSGAYILGRCTVAKKLHDGGLDYFEGYSLENWVCLAYFESKFNPMAIYENTREGY
TGFGFLFQMRGSDWCGDHGRNRCHMSCALLNPNEKTIKCAKTIKKGEGMGAWPTWSRYCQYSDTLARW
LDGCKL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_144634

ORF Size: 438 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_144634.4](#)

RefSeq Size: 886 bp

RefSeq ORF: 441 bp

Locus ID: 131375


UniProt ID: [Q96KX0](#)

Cytogenetics: 3p22.1

Protein Families: Secreted Protein

Gene Summary: Lysozymes (see LYZ; MIM 153450), especially C-type lysozymes, are well-recognized bacteriolytic factors widely distributed in the animal kingdom and play a mainly protective role in host defense. LYZL4 is a member of a family of lysozyme-like genes (Zhang et al., 2005 [PubMed 16014814]).[supplied by OMIM, Apr 2009]

Product images:



Circular map for RG205478