

Product datasheet for **RG205522**

CTRP7 (C1QTNF7) (NM_031911) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTRP7 (C1QTNF7) (NM_031911) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	CTRP7
Synonyms:	CTRP7; ZACRP7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG205522 representing NM_031911 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTTGTCTTGCTCTATGTTACAAGTTTTGCCATTTGTGCCAGTGGACAACCCCGGGTAATCAGTTGA
AAGGAGAGAACTACTCCCCAGGTATATCTGCAGATTCTGGCTTGCCTGGACCTCCAGGGCCCCCTGG
AGCAAATGGTTCCCTGGGCCCATGGTCGCATCGGCCTCCAGGAAGAGATGGTAGAGACGGCAGGAAA
GGAGAGAAAGGTAAAAGGGAAGTGCAGGTTTGTAGAGTAAGACTGGACCGCTAGGTCTTGCCGGTGAGA
AAGGGGACCAAGGAGAGACTGGGAAGAAAGGACCCATAGGACCAGAGGGAGAGAAAGGAGAAGTAGGTCC
AATTGGTCTCTGGACCAAGGGAGACAGAGGAGAACAAGGGGACCCGGGGCTGCCTGGAGTTTGCAGA
TGTGGAAGCATCGTGCTCAAATCCGCCTTTTCTGTTGGCATCACAAACCAGCTACCCAGAAGAAAGACTAC
CTATTATATTTAACAAGTCTCTTCAACGAGGGAGAGCACTACAACCCTGCCACAGGGAAGTTCATCTG
TGCTTTCCAGGGATCTATTACTTTTCTTATGATATCACATTGGCTAATAAGCATCTGGCAATCGGACTG
GTACACAATGGGCAATACCGGATAAAGACCTTCGACGCCAACACAGGAAACCATGATGTGGCTTCGGGGT
CCACAGTCATCTATCTGCAGCCAGAAGATGAAGTCTGGCTGGAGATTTTCTTACAGACCAGAATGGCCT
CTTCTCAGACCCAGGTTGGCAGACAGCTTATTCTCCGGGTTTCTTATACGTTGACACAGATTACCTA
GATTCATATCAGAAGATGATGAATTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG205522 representing NM_031911
Red=Cloning site Green=Tags(s)

MFVLLYVTSFAICASGQPRGNQLKGENYSPRYICSIPLPGPPGPPGANGSPGPHGRIGLPGRDGRDGRK
 GEKGEKGTAGLRGKTGPLGLAGEKGDQGETGKKGPIGPEGEKGEVGPIGPPGPKGDRGEQGDPLPGVCR
 CGSIVLKSAFVSVGITTSYPEERLPPIIFNKVLFNEGEHYNPATGKFCFAFPGIYYFSYDITLANKHLAIGL
 VHNGQYRIKTFDANTGNHDVASGSTVIYLPQPEDEVWLEIFFTDQNGLFSDPGWADSLFSGFLLYVDTDYL
 DSISEDDEL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_031911

ORF Size: 867 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_031911.5](#)

RefSeq Size: 3954 bp

RefSeq ORF: 870 bp

Locus ID: 114905


UniProt ID: [Q9BXJ2](#)

Cytogenetics: 4p15.32

Domains: C1Q, Collagen

Protein Families: Secreted Protein

Product images:



Circular map for RG205522