

Product datasheet for **RG205547**

UBQLNL (NM_145053) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	UBQLNL
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-AC-GFP (PSI00010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG205547 representing NM_145053
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTGGCATGCCATCTCTCGAACATCCAGGATGTCCAGAGTGGATGTCCATCAGGTCTGCTGGCAGACA
AAAATATCTCTTCAAGTGCCACTCGAGTGATAGTGAAGACTGCAGGCAACCAGAAAGACTTTATGGTAGC
TGATGACATCTCGGTAAGGCAGTTCAGGAGATGCTATTGGCTCACTTCCAATGCCAGATGGACCAACTA
GTGCTGGTCTTCATGGGTTGCCTTCTCAAAGACCATGACACACTGAGCCAGAGGGGCATCATGGATGGCC
ACACCATCTACTTGGTCATCAAGTCCAAGCAGGGCTCCAGATCTCTAGCCCATTCCTTCCGGGACCTGCC
AACGAATGATCCCTGCCACCGGGACAGAAACACCAAAGGAAACAGCAGCAGAGTGCACCAACCAACTGGT
ATGAATCAAGCTCCAGTGGAAGTGGCCACTTTGTGGGTCTGATGCACCCAAAGTGCATACCCAAACT
TGGAAGTGAGCCACCCAGAGCGCAAAGCACAGATGCTGGAGAATCTAGCATCCAGCGGCTTCTGTCCAA
CATGGAGTTCATGTGGCAGTTCATTTCAGAACATCTAGACACGCAACAATTGATGCAGCAGAACCCCGAA
GTTTCCCGCTTCTTCTGATAATTCTGAGATCCTATTGCAGACTCTGGAGCTGGCCAGGAACCTTGCTA
TGATCCAAGAGATAATGCAGATCCAACAACCTTCACAAAACCTTGAGTATCCACTGAACCCACAGCCATA
TCTGGGCTTAGAGACAATGCCAGGTGGGAATAATGCCCTGGGTGAGAACCATGCTGATATCAATGATCAA
ATGCTGAACAGCATGCAAGATCCTTTTGAGGAAACCTTTACAGCTCTCCTGGCAGGACAAGTGCTAG
AACAAGTCCAGTCTTACCCCCACCTCCACCACCATCACAGGAACAACAAGACCAGCTCACACAGCATCC
TGCAACCCGAGTCATCTATAATAGCTCTGGTGGTTTCTCTTCAAACACCTCAGCCAATGACACCCCTTAAC
AAGGTCAACCACTTCCAAAGCCAACACTGCTATGATTTCCACCAAGGGCCAGAGCCATATCTGTGCCA
CTCGGCAGCCAGCTGGGATACGAGCTTACCTAGCATAGAGCTTACCCAGCAGCTTCAAGAAGAATACAA
GGATGCCACTGTTTCTAAGTAGCTCCAGACAGACATTAAGGGTGATCTCCAGCTGTCAGATGAACAG
AGCAGCTCCAGATCACAGGAGCATGATGCAGTTGCTTATGAACAACCCCTACCTGGCAGCTCAGATTA
TGTTGTTTCAAGTATGCCCCAGCTGAGTGAACAGTGGAGGCAGCAGCTGCCACATTCCTGCAGCAGAC
ACAGATTTCTGATCTGCTTAGTGCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA




Protein Sequence: >RG205547 representing NM_145053
Red=Cloning site Green=Tags(s)

MWHAISRTSRMSQSGCPSGLLADKNISSSATRVIVKTAGNQKDFMVADDISVRQFKEMLLAHFQCQMDQL
VLVFMGCLLKDHDTLSQRGIMDGHTIYLVIKSQGSRSLAHSFRDLPTNDPCHRDRTKGNSSRVHQPTG
MNQAPVELAHFVGSDAPKVHTQNLEVSHPERKAQMLENPSIQRLLSNMEFMWQFI SEHLDTQQLMQNP
VSRLLLDNSEILLQTLELARNLAMIQEIMQIQPSQNL EYPLNPQPYLGLETMPGGNNALGQNHADINDQ
MLNSMQDPFGGNPFTALLAGQVLEQVQSSPPPPPSQEQQDQLTQHPATRVIIYNSGGFSSNTSANDTLN
KVNHTSKANTAMISTKGQSHICATRQPAIGPALPSIELTQQLQEEYKDATVSLSSSRQTLKGDLQLSDEQ
SSSQITGMMQLLMNPYLAQIMLFTSMPQLSEQWRQQLPTFLQQTQISDLLSA

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_145053

ORF Size: 1425 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_145053.3](#), [NP_659490.3](#)

RefSeq Size: 2319 bp

RefSeq ORF: 1428 bp

Locus ID: 143630


UniProt ID: [Q8IYU4](#)

Cytogenetics: 11p15.4

Domains: UBQ

Protein Families: Druggable Genome

Product images:



Circular map for RG205547