

Product datasheet for **RG205585**

DENND1B (NM_144977) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DENND1B (NM_144977) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	DENND1B
Synonyms:	C1ORF18; C1orf218; FAM31B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG205585 representing NM_144977 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**

ATGGCGGCAGCTCCACGGGAGGAGAAAAGATGGCCCCAACCTGTATTTTCGAATCCTGTGGTATTGTGGA
 AATTCCCAGAGGACTTTGGAGACCAGGAAATACTACAGAGTGTGCCAAAGTTCTGTTTTCCCTTTGACGT
 TGAAAGGGTGTCTCAGAATCAAGTTGGACAGCACTTTACCTTTGTACTGACAGACATTGAAAGTAAACAG
 AGATTTGGATTCTGCAGACTGACGTGAGGAGGCACAATTTGTTTATGCATCCTTAGTTACCTTCCCTGGT
 TTGAAGTGTATTACAAGCTTCTAAATACTCTTGACAGATTACTGGCTAAGGAACTGGAAAATGATTTGAA
 TGAAGTCTCAGATCACTGTATAACCAACCCAGTACCAAAGGCAAATACTCCTGTAATTTGAGTGTGCAT
 TCTACTTCATTGCCCTGATGTAAGTGGACTCCCAACAATACCCGAGAGTAGAAATCTTACAGAATATT
 TTGTTGCCGTGGATGTGAACAACATGCTGCAGCTGTATGCCAGTATGCTGCATGAAAGGCGCATCGTGAT
 TATCTCGAGCAAATTAAGCACTTTAACTGCCTGTATCCATGGATCAGCTGCTCTTCTATACCCAATGTAT
 TGGCAACACATATACATCCAGTGCTTCTCCACACCTGCTGGACTACTGCTGTGCCCAATGCCATACC
 TGATTGGAATACACTCCAGCCTCATAGAGAGAGTGAAAAACAAATCATTGGAAGATGTTGTTATGTTAAA
 TGTTGATACAAACACATTAGAATCACCATTAGTGACTTGAACAACCTACCAAGTATGTGGTCTCGGCC
 TTGAAAAATAAACTGAAGAAGCAGTCTACAGCTACGGGTGATGGAGTAGCTAGGGCCTTTCTTAGAGCAC
 AGGCTGCTTTGTTGGATACTACAGAGATGCACTGAGATACAAACCTGGTGAGCCCATCACTTTCTGTGA
 GGAGAGTTTTGTAAAGCACCGCTCAAGCGTGATGAAACAGTTCCTGGAAGTGCATTAACCTCCAGCTT
 TTTAAGCAGTTTATCGATGGTCGACTGGCAAACTAAATGCAGGAAGGGTTTCTCTGATGTATTGGAAG
 AAGAGTCACTTCAGGTGGCTTTTGTGGAGGTAAAGACAAGTTACAATATGATTATCCATTTTCTCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG205585 representing NM_144977
 Red=Cloning site Green=Tags(s)


MAAAPREEKRWPQPVFSNPVVLWKFPEDFGDQEILQSVPKFCFPFDVERVSQNQVGQHFTFVLTDIESKQ
 RFGFCRLTSGGTICLCILSYLPWFVEVYKLLNTLADYLAKELENDLNETLRSLYNHPVPKANTPVNLSVH
 SYFIAPDVTGLPTIPESRNLTEYFVAVDVNNMLQLYASMLHERRIVIISSKLSTLTACIHGSAALLYPMY
 WQHIYIPVLPPLLDYCCAMPYLIIGIHSSLIERVKNKSLEDVVMNVDNTLTLESPFSDLNNLPSDVVSA
 LKNKLKKQSTATGDGVARAFLRAQAALFGYYRDALRYKPGEPITFCEESFVKHRSSVMKQFLETAINLQL
 FKQFIDGRLAKLNAGRGFSDVFEEIITSGGFCGGKDKLQDYDYPFSQ

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_144977

ORF Size: 1188 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_144977.2](#)

RefSeq Size: 2213 bp

RefSeq ORF: 1281 bp

Locus ID: 163486


UniProt ID: [Q6P3S1](#)

Cytogenetics: 1q31.3

Protein Families: Transmembrane

Gene Summary: Clathrin (see MIM 118955)-mediated endocytosis is a major mechanism for internalization of proteins and lipids. Members of the connectenn family, such as DENND1B, function as guanine nucleotide exchange factors (GEFs) for the early endosomal small GTPase RAB35 (MIM 604199) and bind to clathrin and clathrin adaptor protein-2 (AP2; see MIM 601024). Thus, connectenns link RAB35 activation with the clathrin machinery (Marat and McPherson, 2010 [PubMed 20154091]).[supplied by OMIM, Nov 2010]

Product images:



Circular map for RG205585