

Product datasheet for **RG205696**

TRIM5 alpha (TRIM5) (NM_033092) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM5 alpha (TRIM5) (NM_033092) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	TRIM5 alpha
Synonyms:	RNF88; TRIM5alpha
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG205696 representing NM_033092 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTCTGGAATCCTGGTTAATGTAAGGAGGAGGTGACCTGCCCCATCTGCCTGGAACCTCTGACAC
AACCCCTGAGCCTGGACTGCGGCCACAGCTTCTGCCAAGCATGCCTCACTGCAAACCACAAGAAGTCCAT
GCTAGACAAAGGAGAGAGTAGCTGCCCTGTGTCCGGATCAGTTACCAGCCTGAGAACATACGGCCTAAT
CGGCATGTAGCCAACATAGTGGAGAAGCTCAGGGAGGTCAAGTTAGCCCAGAGGGGCAGAAAGTTGATC
ATTGTGCACGCCATGGAGAGAACTTCTACTCTTCTGTCAGGAGGACGGGAAGGTCATTTGCTGGCTTTG
TGAGCGGTCTCAGGAGCACCGTGGTCACCACACGTTCTCACAGAGGAGGTTGCCCGGGAGTACCAAGTG
AAGCTCCAGGCAGCTCTGGAGATGCTGAGGCAGAAGCAGCAGGAAGCTGAAGAGTTGGAAGCTGACATCA
GAGAAGAGAAAGCTTCTGGAAGACTCAAATACAGTATGACAAAACCAACGTCTTGGCAGATTTTGAGCA
ACTGAGAGACATCCTGGACTGGGAGGAGAGCAATGAGCTGCAAAACCTGGAGAAGGAGGAGGAAGACATT
CTGAAAAGCCTTACGAACTCTGAACTGAGATGGTGCAGCAGACCCAGTCCCTGAGAGAGCTCATCTCAG
ATCTGGAGCATCGGCTGCAGGGGTCAAGTATGGAGCTGCTTCAGGGTGTGGATGGCGTCATAAAAGGAC
GGAGAAGCTGACCTTGAAGAAGCCAGAACTTTTCAAAAAATCAAAGGAGAGTGTTCGAGCTCCTGAT
CTGAAAGGAATGCTAGAAGTGTTCAGAGAGCTGACAGATGTCCGACGCTACTGGGTAAGGAGAAGTCAC
ATTATCATAAGCCACCTGCGGCTTATCATTATTATTATCTTTTATCTTTTAGAATTTTATGTTCTCTATT
AGGCTCATGTTTTAAGATTTATGATTCTCCTTCCAAGACACATAAATTACCCCTCCTTA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG205696 representing NM_033092
 Red=Cloning site Green=Tags(s)

MASGILVNVKEEVTCPICLELLTQPLSLDCGHSFCQACL TANHKKSM LDKGESSCPVCRISYQPENIRPN
 RHVANIVEKLREVKLSPEGQKVDHCAHGEKLLLFQEDGKVICWL CERSQEHARGHHTFL TEEVAREYQV
 KLQAALEMLRQKQQAEELEADIREEKASWKTQIQYDKT NVLADFEQLRDILDWEE SNELQNLEKEEDI
 LKSLTNSETEMVQQTQSLRELISDLEHRLQGSVMELLQGV DGVIKRTENVTLKKPETFPKNQRRVFRAPD
 LKGMLEVFRELTDVRRYWGKEKSHYHKPPCGLSLLL SLSFRILCSLLGSCFKIYDSPSKTHITYPSL

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_033092

ORF Size: 1041 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_033092.1](#), [NP_149083.1](#)

RefSeq Size: 1686 bp

RefSeq ORF: 1044 bp

Locus ID: 85363

UniProt ID: [Q9C035](#)

Cytogenetics: 11p15.4

Domains: zf-B_box, RING

Protein Families: Druggable Genome

Gene Summary: The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein forms homo-oligomers via the coiled-coil region and localizes to cytoplasmic bodies. It appears to function as a E3 ubiquitin-ligase and ubiquitinates itself to regulate its subcellular localization. It may play a role in retroviral restriction. Multiple alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Dec 2009]

Product images:



Circular map for RG205696