

## Product datasheet for **RG205792**

### ABHD12B (NM\_181533) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ABHD12B (NM\_181533) Human Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** ABHD12B  
**Synonyms:** BEM46L3; C14orf29; c14\_5314  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >RG205792 representing NM\_181533  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCTAGGGATCTGGCACACAGTCCCCAGCTGCCGGGGGAAGATGCCAAGGGGAAGGACTGTTGCTGGT  
ATGAAGCAGCCCTTCGTGATGGGAACCAATTATTGTTTATCTTCATGGCAGTGCAGAACACAGGGCAGC  
TTCGCACAGACTGAAGCTGGTAAAGGTGCTGAGTATGGTGGCTTTCATGTCTTGTCTGTTGACTACAGA  
GGATTTGGGACTCTACAGTAAGCCACAGAGGAGGACTGACTACGGATGCCATTTGTGTCTATGAGT  
GGACCAAGGCAAGAAGTGGCATCACTCCCGTGTGTCTCTGGGCCACTCTCTGGGTACAGGAGTTGCAAC  
AAATGCTGCAAAAAGTGTAGAAGAAAAAGGATGCCAGTTGATGCTATTGTCTTGAAGCTCCATTTACC  
AACATGTGGGTTGCAAGTATCAATTATCCCTTGTAAAGATTTACCGAACATTCCAGGATTTTTACGTA  
CACTTATGGATGCCCTGAGAAAAGACAAAATAATCTTTCCTAATGATGAAAATGTTAAATTCCTTCTTC  
TCCTCTTCTCATCTTACATGGAGAGGATGACAGGACAGTGCCTTTGGAGTATGGGAAAAGCTCTATGAA  
ATTGCACGCAATGCATACAGGAACAAAGAGAGGGTCAAGATGGTTATCTTTCCTCTGGCTTCCAACACA  
ACCTGCTTGTAAAAGCCCCACACTGTTAATAACCGTGAGAGATTTCTGAGCAAGCAGTGGTCA

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >RG205792 representing NM\_181533  
Red=Cloning site Green=Tags(s)

MLGIWHTVPSCRGEDAKGKDCCWYEAALRDGNPIIVYLHGSAEHRAASHRLKLVKVLSDGGFHVLSVDYR  
 GFGDSTGKPTTEGLTTDAICVYEWTKARSGITPVCLWGHSLGTGVATNAAKVLEEKGCPVDAIVLEAPFT  
 NMWVASINYPLLKLYRNIPGFLRTLMDALRKDKIIFPNDENVKFLSSPLLILHGEDDRTVPLEYGKKLYE  
 IARNAYRNKERVKMVIFFPPGFQHNLLCKSPTLLITVRDFLSKQWS

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_181533

**ORF Size:** 765 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_181533.3](#), [NP\\_853511.2](#)

**RefSeq Size:** 1580 bp


**RefSeq ORF:** 768 bp

**Locus ID:** 145447

**UniProt ID:** [Q7Z5M8](#)

**Cytogenetics:** 14q22.1

**Product images:**



Circular map for RG205792