

Product datasheet for **RG205822**

PGAP6 (NM_021259) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PGAP6 (NM_021259) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	PGAP6
Synonyms:	GPI-PLA2; M83; TMEM6; TMEM8; TMEM8A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RG205822 representing NM_021259
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**C

ATGGGCCGGGCTGGCACCAGGACCGGGGGCGAGCGGTGGCCGCGGTGGTGGCGGGCCGCTGCTGCTGC
 TGTGCTTGCCCCGGCCCCGCTGCCTCCGCGGCTACAGCGGAAGAGCGAGGTGGGGCTGGTGTCCGA
 GCACTTCTCACAGGCCCCGAGAGGCTGTCTTCTACAGCTGGTATGGCAGTGCCAGGCTCTTCCGCTTC
 CGCGTGCCCCAGATGCTGTGCTTCTACGCTGGCTCCTGCAGGTCTCCCGGAGAGCGGCGCTGCCTGCA
 CCGACGCGGAGATCACCGTGCACCTCCGTTCCGGCGCCCTCCGGTCATCAACCCGCTGGGACACAGCTT
 CCCGGACGACACCGCGGTGCAGCCCTCCTTCCAGGTGGGGTGGCGCTGAGCACCACCGAGAAGCAAT
 GCCTCCGTCAACGTTTCCACCCGGCCCCGGGACTGGTTCGTGGCGCCACCTGCCCCCTCATCCC
 AGAAGATCGAGTTGAAGGGCTTGGCTCCACCTGTGCCTACGTCTTCCAGCCTGAAGTGTGGTCACGCG
 GGTGGTCGAGATTCCATCATGGAGCCGACGTGCCCCCTCCTCAGACCCTCCTCTCCATCCCAGCTAC
 CTCAAGGTCTTTGTCCCGATTACACGCGGGAGCTTCTGCTGGAGCTGCGGGACTGCGTGTCCAATGGGA
 GCCTGGGCTGCCCGTGCGTCTACCGTGGGCCCGGTACCCCTGCCTAGCAACTTCCAGAAGGTGCTCAC
 CTGCACCGGTGCCCTGGCCCTGCCGCTGCTGCTGCCCTACCGCCCTGGGACCGGTGGTGTCAAGTG
 ACAGCTGAGAGCTGGTGGGGCCCTCGGGACAGTGGCTTTCAGTGCTGTAGTGCCCTCACAGCTTGCA
 GGCCACGAGCGGTGACCGTCCAGCCCCCTTCTGCAGAGCAGCAAAACAGAGCTTCAATGCCTCCTCTGG
 TCTGCTGTCCCCGAGCCCCGACCACAGGACCTGGGACAGAGTGGCAGGGTGGACCGCAGCCCCCTTCTGC
 CTCACAACTACCCAGTCACGCGGGAGGACATGGACGTGGTGTGGTGCACCTTCCAGCCCTGGACAGGG
 TCTCGGTGAGGGTGTGTTCCGACACGCCCTCCGTGATGCGGCTGCGCCTGAACACCGCATGGACAGCGG
 GGGTTCCTCACCATCTCCCTGCGGGCCAACAAGACAGAGATGCGGAACGAGACCGTCGTAGTGGCCTGC
 GTGAATGCTGCCTCGCCCTTCTTGGCTTCAATACTTCGCTCAACTGCACCACAGCCTTCTTCCAGGGCT
 ACCCTTTGTCTCTGAGCGCTGGTCTCGCAGGGCCAACCTCATCATCCCCTACCCAGAGACAGAACTG
 GTACCTCTCCCTGCAGCTCATGTGCCCTGAGAATGCTGAGGACTGTGAGCAGGCTGTGGTCCACGTGGAG
 ACCACCTTGACCTGGTGCCTGTTTGAACGATTGTGGACCTATGGCCAGTGCCCTCCTGCTCCGACAGC
 ACAGCTACCTGTATGCCAGCTGCAGCTGCAAGGCAGGCTGGCGTGGTGGAGCTGCACGGACAACAGCAC
 AGCCACAGCGGTGGCCAGCAGAGGGCGGCCACACTGCTGCTCACGCTCAGCAACCTCATGTTCTTGCC
 CCCATCGCCGTCTCAGTGCAGGCTTCTTCTGTTGGAGGCTCCGTCTACGCTACACCATGTTCTTCT
 CCACGTTCTACCACGCTGCGACCGAGCCCGGGAGGCGGTGCTGTGCATCCTCAGCTACGACACGCTGCA
 GTACTGCGACTTCTTGGGCTCCGGGCGGCCATCTGGGTACCATCCTGTGCATGGCAGCGCTCAAGACA
 GTCCTGAAATACGTGCTGTTTCTTCTGGGTACACTGGTATCGCCATGTCCTTGACGCTGGACCGCAGGG
 GCATGTGGAACATGCTGGGGCCCTGCCTCTTGCCTTCGTGATCATGGCTCCATGTGGGCTTACCGCTG
 CGGGACCGGCGCCAGTGCTACCCACCTCGTGGCAGCGCTGGGCTTCTACCTCCTGCCCGGCTCTCT
 ATGGCCTCTGTGGGCATCGCCATCTACACCTCCATGATGACTAGCGCAACTACTACTACCCACAGCA
 TCTGGCACATCCTGCTGGCCGGGAGCGCAGCCTTGCTGCTGCCGCCACCTGACGAGCCCGCGAGCCCTG
 GGCCTGCTCGCAGAAATCCCCTGCCACTATCAGATCTGCAAGAACGATCGGGAGGAAGTGTACGCAAGT
 ACG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA


Protein Sequence: >RG205822 representing NM_021259
 Red=Cloning site Green=Tags(s)

MGRAGTGTGGEAAVAVAGPLLLLLLARPPASAGYSGKSEVGLVSEHFSQAPQRLSFYSWYGSARLFRF
 RVPPDAVLLRWLLQVSRESGAACDAEITVHFRSGAPPVINPLGTSFPDDTAVQPSFQVGVPLSTAPRSN
 ASVNVSHAPGDWFAAHLPPSSQKIELKGLAPTCAYVFQPELLVTRVVEISIMEPDVPLPQTLLSHPSY
 LKVFVPDYTRELLLELRDCVSNGLGCPVRLTVGPVTLPSNFQKVLCTGAPWPCRLLLPSPWDRWLQV
 TAESLVGPLGTVAFSAVAALTACRPRSVTVQPLLQSSQNSFNASSGLLSPSPDHQDLGRSGRVDSPFC
 LTNYPVTREMDVVSVHFQPLDRVSVRVCSDTPSVMRLRLNTGMDSGGSLTISLRANKTEMRNETVVVAC
 VNAASPFLGFNTSLNCTTAFQGYPLSLAWSRRANLIIPYPETDNWYLSQLMCPENAEDCEQAVVHVE
 TTLYLVPCLNDCGPYGQCLLLRRHSYLYASCCKAGWRGWSCTDNSTAQTVAQQRAATLLTLNLMFLA
 PIAVSVRRFFLVEASVYAYTMFFSTFYHACDQPGEA VLCILSYDTLQYCDFLGSGAAIWVILCMARLKT
 VLKYVLFLLGTLVIAMSLQLDRGMWMLGPCLFAFVIMASMWAYRCGHRRCYPTSWQRWAFYLLPGVS
 MASVGIAIYTSMMTSDNYYYTHSIWHILLAGSAALLPPPDQPAEPWACSQKFPCHYQICKNDREEL YAV
 T

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_021259


ORF Size: 2313 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_021259.1, NP_067082.1</u>
RefSeq Size:	2561 bp
RefSeq ORF:	2316 bp
Locus ID:	58986
UniProt ID:	<u>Q9HCN3</u>
Cytogenetics:	16p13.3
Protein Families:	Transmembrane
Gene Summary:	Involved in the lipid remodeling steps of GPI-anchor maturation. Lipid remodeling steps consist in the generation of 2 saturated fatty chains at the sn-2 position of GPI-anchor proteins (GPI-AP). Has phospholipase A2 activity that removes an acyl-chain at the sn-2 position of GPI-anchors during the remodeling of GPI. Required for the shedding of the GPI-AP TDGF1, but not CFC1, at the cell surface. Shedding of TDGF1 modulates Nodal signaling by allowing soluble TDGF1 to act as a Nodal coreceptor on other cells (PubMed:27881714). Also indirectly involved in the translocation of RAC1 from the cytosol to the plasma membrane by maintaining the steady state amount of CAV1-enriched plasma membrane subdomains, stabilizing RAC1 at the plasma membrane (PubMed:27835684). In contrast to myomaker (TMEM8C), has no fusogenic activity (PubMed:26858401).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RG205822