

Product datasheet for **RG205859**

GRP78 (HSPA5) (NM_005347) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GRP78 (HSPA5) (NM_005347) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	GRP78
Synonyms:	BIP; GRP78; HEL-S-89n
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RG205859 representing NM_005347
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGCATCGCC

ATGAAGCTCTCCCTGGTGCCGCGATGCTGCTGCTCAGCGCGCCGGGCCGAGGAGGACAAGA
AGGAGGACGTGGGCACGGTGGTCGCGATCGACCTGGGGACCACCTACTCCTGCGTCGGCGTGTCAAGAA
CGGCCGCGTGGAGATCATCGCCAACGATCAGGGCAACCGCATCACGCCGTCCTATGTCGCCCTTCACTCCT
GAAGGGGAACGTCTGATTGGCGATGCCGCCAAGAACCAGCTCACCTCCAACCCGAGAACACGGTCTTTG
ACGCCAAGCGGCTCATCGGCCGACGTGGAATGACCCGCTGTGCAGCAGGACATCAAGTTCTTGCCGTT
CAAGGTGGTTGAAAAGAAAATAAACCATACATTCAGTTGATATTGGAGGTGGGCAAACAAGACATTT
GCTCCTGAAGAAATTTCTGCCATGGTTCTCACTAAAATGAAAGAAACCGCTGAGGCTTATTTGGAAAGA
AGGTTACCCATGCAGTTGTTACTGTACCAGCCTATTTAATGATGCCAACGCCAAGCAACCAAAGACGC
TGGAACTATTGCTGGCCTAAATGTTATGAGGATCATCAACGAGCCTACGGCAGCTGCTATTGCTTATGGC
CTGGATAAGAGGGAGGGGAGAAGAATCCTGGTGTGGACCTGGGTGGCGGAACCTTCGATGTGTCTC
TTCTCACCATTGACAATGGTGTCTTCAAGTTGTGGCCACTAATGGAGATACTCATCTGGGTGGAGAAGA
CTTTGACCAGCGTGTGATGGAACACTTCATCAAACTGTACAAAAAGAACGGGCAAAGATGTCAGGAAA
GACAATAGAGCTGTGCAGAACTCCGGCGCAGGTAGAAAAGGCCAAACGGGCCCTGTCTTCTCAGCATC
AAGCAAGAATTGAAATGAGTCTTCTATGAAGGAGAAGACTTTTCTGAGACCCTGACTCGGGCCAAATT
TGAAGAGCTCAACATGGATCTGTCCGGTCTACTATGAAGCCCGTCCAGAAAGTGTGGAAGATTCTGAT
TTGAAGAAGTCTGATATTGATGAAATGTTCTTGTGGTGGCTCGACTCGAATTCCAAAGATTGAGCAAC
TGGTAAAGAGTTCTTCAATGGCAAGGAACCATCCCGTGGCATAAACCAGATGAAGCTGTAGCGTATGG
TGCTGCTGTCCAGGCTGGTGTGCTCTCTGGTGATCAAGATACAGGTGACCTGGTACTGCTTGTATGT
CCCCTTACACTTGGTATTGAAACTGTGGGAGGTGTGATGACCAAACCTGATTCCAAGGAACACAGTGGTGC
CTACCAAGAAGTCTCAGATCTTTTCTACAGCTTCTGATAATCAACCAACTGTTACAATCAAGGTCTATGA
AGGTGAAAGACCCTGACAAAAGACAATCATCTTCTGGGTACATTTGATCTGACTGGAATTCCTCCTGCT
CCTCGTGGGTCCACAGATTGAAGTCACTTTGAGATAGATGTGAATGGTATTCTTCGAGTGACAGCTG
AAGACAAGGTACAGGGAACAAAATAAGATCACAATCACAATGACCAGAATCGCCTGACACCTGAAGA
AATCGAAAGGATGGTTAATGATGCTGAGAAGTTGCTGAGGAAGACAAAAGCTCAAGGAGCGCATTGAT
ACTAGAAATGAGTTGAAAGCTATGCCTATTCTCTAAAGAAATCAGATTGGAGATAAAGAAAAGCTGGGAG
GTAACCTTCTCTGAAGATAAGGAGACCATGGAAAAAGCTGTAGAAGAAAAGATTGAATGGCTGGAAAG
CCACCAAGATGCTGACATTGAAGACTTCAAGCTAAGAAGAAGGAACTGGAAAGAAATTGTTCAACCAATT
ATCAGCAAACCTATGGAAGTGCAGGCCCTCCCCAACTGGTGAAGAGGATACAGCAGAAAAAGATGAGT
TG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG205859 representing NM_005347
Red=Cloning site Green=Tags(s)

```

MKLSLVAAMLLLLSAARAEEDKKEDVGTVVGIDLGTTYSCVGVFKNGRVEIIANDQGNRITPSYVAFTP
EGERLIGDAAKNQLTSNPENTVFDKRLIGRTWNDPSVQQDIKFLPFKVVKKTKPYIQVDIGGGQTKTF
APEEISAMVLTKMKETAAYLGGKVTTHAVVTPAYFNDAQRQATKDAGTIAGLNVMRIINEPTAAAIAYG
LDKREGEKNILVFDLGGGTFDVSLLTIDNGVFEVVATNGDTHLGGEDFDQRMVMEHFIKLYKKKTGKDVRK
DNRAVQKLRREVEKAKRALSSQHARIEIESFYEGEDFSETLTRAKFEELNMDLFRSTMKPVQKVLSD
LKKSDIDEIVLVGGSTRIPKIQQLVKEFFNGKEPSRGINPDEAVAYGAAVQAGVLSGDQDTGDLVLLDVC
PLTLGIETVGGVMTKLIPRNTVVPTKKSQIFSTASDNQPTVTIKVYEGERPLTKDNHLLGTFDLTGIPPA
PRGVPQIEVTFEIDVNGILRVTAEDKGTGNKNKIIITNDQNRLTPEEIERMVNDAEKFAEEDKKKERID
TRNELESYAYSLKNQIGDKEKLGKLSSEDKETMEKAVEEKIEWLESHQDADIEDFKAKKKELEEIVQPI
ISKLYGSAGPPPTGEEDTAEKDEL
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_005347

ORF Size: 1962 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005347.5](#)

RefSeq Size: 3925 bp

RefSeq ORF: 1965 bp

Locus ID: 3309

UniProt ID: [P11021](#)

Cytogenetics: 9q33.3


Domains: HSP70

Protein Families: Druggable Genome

Protein Pathways: Antigen processing and presentation, Prion diseases

Gene Summary: The protein encoded by this gene is a member of the heat shock protein 70 (HSP70) family. This protein localizes to the lumen of the endoplasmic reticulum (ER) where it operates as a typical HSP70 chaperone involved in the folding and assembly of proteins in the ER and is a master regulator of ER homeostasis. During cellular stress, as during viral infection or tumorigenesis, this protein interacts with the transmembrane stress sensor proteins PERK (protein kinase R-like endoplasmic reticulum kinase), IRE1 (inositol-requiring kinase 1), and ATF6 (activating transcription factor 6) where it acts as a repressor of the unfolded protein response (UPR) and also plays a role in cellular apoptosis and senescence. Elevated expression and atypical translocation of this protein to the cell surface has been reported in viral infections and some types of cancer cells. At the cell surface this protein may facilitate viral attachment and entry to host cells. This gene is a therapeutic target for the treatment of coronavirus diseases and chemoresistant cancers. [provided by RefSeq, Jul 2020]

Product images:



Circular map for RG205859