

Product datasheet for **RG206022**

Kallikrein 11 (KLK11) (NM_006853) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kallikrein 11 (KLK11) (NM_006853) Human Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kallikrein 11
Synonyms:	PRSS20; TLSP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>RG206022 representing NM_006853 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGATTCTGCAGTTAATCCTGCTTGTCTGGCAACAGGGCTTGTAGGGGGAGAGACCAGGATCATCA
AGGGGTTTCGAGTGCAAGCCTCACTCCCAGCCCTGGCAGGCAGCCCTGTTTCGAGAAGACGGGCTACTCTG
TGGGGCGAGCTCATCGCCCCAGATGGCTCTGACAGCAGCCCACTGCCTCAAGCCCCGTACATAGTT
CACCTGGGGCAGCACAACTCCAGAAGGAGGAGGGCTGTGAGCAGACCCGGACAGCCACTGAGTCCCTCC
CCCACCCCGGCTTCAACAACAGCCTCCCAACAAGACCACCGCAATGACATCATGCTGGTGAAGATGGC
ATCGCCAGTCTCCATCACCTGGGCTGTGCGACCCCTCACCTCTCCTCACGCTGTGTCACTGCTGGCACC
AGCTGCCTCATTTCCGGCTGGGGCAGCAGTCCAGCCCCAGTTACGCCTGCCTCACACCTTGCATGCG
CCAACATCACCATCATTGAGCACCAGAAGTGTGAGAACGCCTACCCCGGCAACATCACAGACCCATGGT
GTGTGCCAGCGTGCAGGAAGGGGCAAGGACTCCTGCCAGGGTACTCCGGGGCCCTCTGGTCTGTAAC
CAGTCTCTTCAAGGCATTATCTCCTGGGGCCAGGATCCGTGTGCGATCACCCGAAAGCCTGGTGTCTACA
CGAAAGTCTGCAAATATGTGGACTGGATCCAGGAGACGATGAAGAACAAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >RG206022 representing NM_006853
Red=Cloning site Green=Tags(s)



MRILQLILLALATGLVGGETRIIKGFECKPHSQPWQAALFEKTRLLCGATLIAPRWLLTAHCLKPRYIV
 HLGQHNLQKEEGCEQTRTATESFPHPGFNNSLPNKDHRNDIMLVKMASPVSITWAVRPLTLSSRCVTAGT
 SCLISGWGSTSSPQLRLPHTLRCANITIIIEHQKCNAYPGNITDTMVCASVQEGGKDCSQGDSGGPLVCN
 QSLQGIISWGDPCAITRKPGVYTKVCKYVDWIQETMKNN

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_006853

ORF Size: 750 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_006853.3](#)

RefSeq Size: 1207 bp

RefSeq ORF: 753 bp

Locus ID: 11012

UniProt ID: [Q9UBX7](#)


Cytogenetics: 19q13.41

Domains: Tryp_SPc

Protein Families: Druggable Genome, Protease, Secreted Protein

Gene Summary: Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Alternate splicing and the use of alternate promoters results in multiple transcript variants encoding distinct isoforms which are differentially expressed. [provided by RefSeq, Dec 2016]

Product images:



Circular map for RG206022