

Product datasheet for **RG206144**

C14orf48 (LINC00521) (NM_152777) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: C14orf48 (LINC00521) (NM_152777) Human Tagged ORF Clone

Tag: TurboGFP

Symbol: C14orf48

Synonyms: MGC33356, c14_5713

Mammalian Cell Selection: Neomycin

Vector: pCMV6-AC-GFP (PS100010)

E. coli Selection: Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >RG206144 representing NM_152777
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACACAGGCCAGAGAGCTGACCCAAGCAATCCTGGTGACAAGGAAGGGGACCTTCAAGGGCTGTGGC
AGGAATCTACCAGCTCCAGGCTAAGCAGAAGAAGCTCAAGAGAGAAGTCGAGAAGCACAAAGCTTTTGA
AGACTATCTGATTAAGGTCCTTGAGAAAATCCCGAGGGCTGCACGGGATGGGAGGAGCCGGAGGAGGTG
CTGGTGGAGGCCACGGTGAAGCACTACGGGAAGCTCTTACAGCCAGCCAGGACACGCAGAAGCGCCTCG
AGGCCTTCTGCCAGATGATCCAGGCTGTCCACCGAGCCTGGAGTCTCTGGAGGAGACCACAGGGCTCT
CATAGCCTCAAGATCCGGCTGTGTGTCAGTGCAGAAGAAGTGTACCGCAAGCAGGAGCAGTGGTGGCAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >RG206144 representing NM_152777
Red=Cloning site Green=Tags(s)


MDTGQRADPSNPGDKEGDLQGLWQELYQLQAKQKCLKREVEKHKLFDYLIKVLEKIPEGCTGWEEPEEV
LVEATVKHYGKLFTASQDTQKRLEAFCQMIQAVHRSLESLEEDHRALIASRSGCVSCRRSATASRSSGGS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_152777

ORF Size: 420 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_152777.2](#), [NP_689990.1](#)

RefSeq Size: 1290 bp

RefSeq ORF: 422 bp
 Locus ID: 256369
 Cytogenetics: 14q32.12

Product images:



Circular map for RG206144